

## **Air Back**

*Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market.*

The Prism™ Air Back was designed for users who require moderate to moderately complex postural support and offers the most aggressive lateral support in the Prism line of products.

Extras such as the choice of bladder position, the ease of attachment and the unique method of adjusting for a wide range of angles, makes this Back the most versatile and truly modular air model on the market.

Combine this with Prism's sleek and attractive tailored look for a Back that will appeal to a wide variety of users!





### SPOTLIGHT ON

- Lightweight, sturdy aluminum shell
- Seat angle adjustment from 82 degrees to 114 degrees optimizing users' positioning
- Attractive cover design and handy pocket to conceal inflation valves
- Cane clamps included
- Adjustable lumbar and lateral bladders for accurate, modifiable support
- Quick push button release for easy installation

### WARRANTY


Three (3) years limited on Back and 90 days on cover, 1 year on hardware.

### PRODUCT WEIGHT

5.3lbs Basic

### USER WEIGHT

250 lbs.



**Stock Sizes—Same day shipments**

16" x 16"	AB1616
18" x 16"	AB1816
20" x 16"	AB2016
16" x 20"	AB1620
18" x 20"	AB1820
20" x 20"	AB2020

**Non-Stock Sizes—Shipments within 3 business days**

14" x 14"	AB1414
17" x 16"	AB1716
19" x 16"	AB1916
16" x 18"	AB1618
18" x 18"	AB1818
17" x 20"	AB1720
19" x 20"	AB1920

*Custom orders will rarely be turned away.....we can accommodate as low as 11" in back height and ship it within 7 days!*



3223 Orlando Drive, Mississauga, ON, L4V 1C5, 905-671-1661, 888-737-4011, fax 866-670-3377 www.future-mobility.com  
 Future Mobility Healthcare Inc. Orion and Prism are trademarks of Future Mobility Healthcare Inc.