

Lunar Back

Future Mobility Healthcare manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Renown for quality and service we bring innovation to all of our products for the home healthcare market.

The Prism™ Lunar Back was designed for users who require aggressive lateral support and more complex positioning.

Extras such as fixed lateral supports or adjustable swing away lateral supports make this back the most versatile and truly modular model on the market.

Combine this with Prism's sleek and attractive tailored look for a back that will appeal to a wide variety of users!



MADE IN CANADA



SPOTLIGHT ON

- Lightweight, sturdy aluminum shell
- Attractive cover design
- Replaceable foam
- Optional Lateral supports adjustable horizontally
- Pins adjustable to accommodate up to 2” difference in seat width
- Available in 3 laterals styles (no lateral, basic lateral, ultra lateral)
- Available in 3 back styles (low back, mid back and full back)
- Optional Lateral supports available in 4 styles (Fixed, Fixed Height Adjustable, Fixed Offset and Adjustable Swingaway)
- Uses newly design Lunar hardware for clamping to tubes from 7/8” diameter to 1-1/8” using no additional inserts

WARRANTY

Three (3) years limited on Back and 90 days on cover, 1 year on hardware.

PRODUCT WEIGHT

4.3 lbs (without lateral supports)
4.9 lbs (with lateral supports)

USER WEIGHT

250 lbs

• *Stock sizes—same day shipments*

LOW BACK STYLE—NO LATERAL

LBSNL1610 16” x 10” LBSNL1810 18” x 10”
LBSNL1612 16” x 12” LBSNL1812 18” x 12” LBSNL2012 20” x 12”

LOW BACK STYLE—BASIC LATERAL

LBSBL1610 16” x 10” LBSBL1810 18” x 10”
LBSBL1612 16” x 12” LBSBL1812 18” x 12” LBSBL2012 20” x 12”

MID STYLE BACK—NO LATERAL

MBSNL1614 16” x 14” MBSNL1814 18” x 14” MBSNL2014 20” x 14”
MBSNL1616 16” x 16” MBSNL1816 18” x 16” MBSNL2016 20” x 16”
MBSNL1618 16” x 18” MBSNL1818 18” x 18” MBSNL2018 20” x 18”

MID STYLE BACK—BASIC LATERAL

MBSBL1614 16” x 14” MBSBL1814 18” x 14” MBSBL2014 20” x 14”
MBSBL1616 16” x 16” MBSBL1816 18” x 16” MBSBL2016 20” x 16”
MBSBL1618 16” x 18” MBSBL1818 18” x 18” MBSBL2018 20” x 18”

MID STYLE BACK—ULTRA LATERAL

MBSUL1614 16” x 14” MBSUL1814 18” x 14” MBSUL2014 20” x 14”
MBSUL1616 16” x 16” MBSUL1816 18” x 16” MBSUL2016 20” x 16”
MBSUL1618 16” x 18” MBSUL1818 18” x 18” MBSUL2018 20” x 18”

FULL STYLE BACK—NO LATERAL

FBSNL1620 16” x 20” FBSNL1820 18” x 20” FBSNL2020 20” x 20”
FBSNL1622 16” x 22” FBSNL1822 18” x 22” FBSNL2022 20” x 22”
FBSNL1824 18” x 24” FBSNL2024 20” x 24”

FULL STYLE BACK—BASIC LATERAL

FBSBL1620 16” x 20” FBSBL1820 18” x 20” FBSBL2020 20” x 20”
FBSBL1622 16” x 22” FBSBL1822 18” x 22” FBSBL2022 20” x 22”
FBSBL1824 18” x 24” FBSBL2024 20” x 24”

FULL STYLE BACK—ULTRA LATERAL

FBSUL1620 16” x 20” FBSUL1820 18” x 20” FBSUL2020 20” x 20”
FBSUL1622 16” x 22” FBSUL1822 18” x 22” FBSUL2022 20” x 22”
FBSUL1824 18” x 24” FBSUL2024 20” x 24”

• *Non-stock sizes—shipments within 3 business days*

LOW BACK STYLE

NO LATERAL LBSNL1410 14” x 10” LBSNL1412 14” x 12”

BASIC LATERAL LBSBL1410 14” x 10” LBSBL1412 14” x 12”

MID STYLE BACK—NO LATERAL

FBSNL1414 14” x 14” FBSNL1416 14” x 16” FBSNL1418 14” x 18”

MID STYLE BACK—BASIC LATERAL

MBSBL1414 14” x 14” MBSBL1416 14” x 16” MBSBL1418 14” x 18”

MID STYLE BACK—ULTRA LATERAL

MBSUL1414 14” x 14” MBSUL1416 14” x 16” MBSUL1418 14” x 18”

FULL STYLE BACK

NO LATERAL FBSNL2222 22” x 22” FBSNL2224 22” x 24”

BASIC LATERAL FBSBL2222 22” x 22” FBSBL2224 22” x 24”

ULTRA LATERAL FBSUL2222 22” x 22” FBSUL2224 22” x 24”

• *Custom sizes—shipments within 7 business days*

LUNARCUS ____W x ____H LATERAL STYLE: NL BL UL

BACK STYLE: LB MB FB

