

New



Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market.



**Attendant or Occupant
Propelled!**

The Aqua Commode is the most innovative, affordable shower commode on the market today. Developed for those individuals requiring assistance for shower or toilet use, the Aqua Commode allows the user to function independently and become more involved in their lifestyles.

The Aqua Commode is a complete over the toilet commode and can also be used near bed side. It has standard back height adjustment, incredible range of seat-to-floor heights and creative streamlined design. It is available in assist and self propelled style.

With its progressive design, the Aqua Commode is sure to become the most popular shower assistance commode on the market.

Future M  bility Healthcare



SPOTLIGHT ON:

- Sleek, cozy and modern in design
- Frame available in steel or stainless steel
- Comfortable, Quick flip back and removable arm rest
- Standard removable swing away foot rest
- Comfortable seating with height adjustable caster stem
- Incredible low 16" seat-to-floor height
- Unique release clip on casters and back for adjustment of seat to floor and back height
- Easy to remove seat and pan
- Height adjustable back
- Available in self propelled style (Optional)
- Spring return push-to-lock wheel lock (Optional)



FEATURES:

- Seat width: 17"
- Seat depth: 17"
- Seat-to-floor heights: 16" - 19"
- Front rigging: 50 degree standard
- 5" casters
- Adjustable foot rest
- 20", 22", 24" rear wheels (Optional)

WARRANTY

Commode frame made from steel is warranted for two(2) years from the original purchase date while stainless steel frame is warranted for five (5) years from the original purchase date. All other parts and hardware are warranted for one(1) year.

SPECIFICATIONS

Chair width: 24"
Chair length: 32.5"
Chair height: 41"
Chair weight: 38 lbs
Capacity : 250 lbs

Future Mobility Healthcare