


Elite Cushion



The Prism Elite Cushion was designed for users that require moderate postural support and offers gentle contours helping to promote correct seating posture.

A top layer of soft, comfortable poly foam, along with the visco and high density polyurethane foam offers pressure distribution and relief. We truly have designed an Elite cushion.



Prism™ Elite Cushion

Stock Sizes—Next Day Shipments

16" x 16"	EC1616	16" x 18"	EC1618
18" x 16"	EC1816	18" x 18"	EC1818
20" x 16"	EC2016	20" x 18"	EC2018

Non Stock Sizes—Shipments within 3 business days

16" x 20"	EC1620	17" x 16"	EC1716
17" x 18"	EC1718	17" x 20"	EC1720
18" x 20"	EC1820	20" x 20"	EC2020

CUSHION THICKNESS

- A—Height at adductor = 4"
- B—Height at abductor = 4 1/2"
- C—Height at rear IT = 3"
- D—Leg trough = 3"

WARRANTY

Two (2) years limited on Cushion and 90 days on cover.

PRODUCT WEIGHT

2.5 lbs.

USER WEIGHT

250 lbs.

SPOTLIGHT ON

- Waterfall front edge for foot propelling
- Lightweight cushion
- Non-slip bottom on washable incontinent cover
- Core layer of visco-elastic foam to provide comfort and pressure reduction

