

Future Mobility Healthcare

Orion Explorer Magellan



Order Form/Price List

(August 6, 2010)

Account Information:

Account #: _____
 Name: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____
 Comment: _____

Date: _____
 PO# _____
 Ship To: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____

FRAME TYPE

Explorer Magellan
 (includes full suspension base, std joystick and lap belt)

EXPLORER \$6,045

FENDER COLOR

Black _____
 Blue _____
 Red _____
 Green _____

FBK No Charge
 FBU No Charge
 FRD No Charge
 FGR No Charge

CONTROLLER

Type

VR2 _____
 RNET _____

VR2 Standard
 RNET \$1,345

Controller Side

Left _____
 Right _____

CLFT No Charge
 CRHT No Charge

Mounting Options

Fixed _____
 Swingaway _____

MOF Standard
 MSWG \$284

Joystick Options

VR2 Joystick _____
 Straight Extension _____
 T Extension _____

JVR2 Standard
 JSE \$58
 JST \$83

DRIVE WHEELS

Pneumatic _____
 Airless _____

PNEU Standard
 ALES \$145

BATTERIES

8G22NF Batteries (pair) _____

8G22NF \$500

CHARGERS

1 4A On Board _____
 8A Off Board _____

CONB Standard
 COFB \$329

STANDARD SEATING

(includes seat frame, solid seat pan, 2-point adjustable height flipback removeable arms)

SEAT WIDTH

2 16" Seat Width _____
 17" Seat Width _____
 18" Seat Width _____
 19" Seat Width _____
 3 20" Seat Width _____
 4 21" Seat Width _____
 4 22" Seat Width _____

STDSEAT \$842

MSW16 No Charge
 MSW17 No Charge
 MSW18 No Charge
 MSW19 No Charge
 MSW20 No Charge
 MSW21 \$274
 MSW22 \$274

Standard Features:

- 1 Full Suspension Base
- 1 Standard Joystick
- 1 Lap Belt
- 2 Pneumatic Wheels
- 2 Battery Box Covers
- 1 Battery Charger
- 1 VR2 Controller
- 2 Fenders
- 4 Caster Wheels
- 2 Motors
- 2 Oval Arm Pads
- 1 One Piece Footboard or
 2 Pin Style Leg Rests
- 1 PRISM Ideal Cushion
- Angle Adjustable Back
 Cane (87° TO 110°)
- Accepts Group 22 Batteries

FOOTNOTES:

1. 8A offboard battery charger provides a faster recharge if overnight is not acceptable
2. seat pan is adjustable from 16" wide to 20" wide
3. seat pan is adjustable from 20" wide to 24" wide
4. upcharge for seat widths 20" - 24" \$274

4	<input type="checkbox"/> 23" Seat Width -----	MSW23	\$274
4	<input type="checkbox"/> 24" Seat Width -----	MSW24	\$274

SEAT DEPTH

5	<input type="checkbox"/> 17" Seat Depth -----	MSD17	No Charge
	<input type="checkbox"/> 18" Seat Depth -----	MSD18	No Charge
	<input type="checkbox"/> 19" Seat Depth -----	MSD19	No Charge
	<input type="checkbox"/> 20" Seat Depth -----	MSD20	No Charge

SEAT-TO-FLOOR HEIGHT

	<input type="checkbox"/> 17" STF-----	SSTF17	No Charge
	<input type="checkbox"/> 18" STF-----	SSTF18	No Charge
	<input type="checkbox"/> 19" STF-----	SSTF19	No Charge

CANE OPTIONS

6.7	<input type="checkbox"/> Cane Back w/handles-----	CBH	Standard
	<input type="checkbox"/> Cane Back no handles-----	CBNH	No Charge

BACK OPTIONS

	<input type="checkbox"/> Orion Back (16", 18", 20") available height only-	OBK	Standard
	<input type="checkbox"/> Prism Basic Back on chair Stock -----	BBOC	\$150
	<input type="checkbox"/> Prism Basic Back on chair Non Stock -----	BBOC	\$175
	<input type="checkbox"/> Prism Basic Back on chair Bariatric -----	BBOC	\$250
	<input type="checkbox"/> Prism Basic Back Custom on chair-----	BBCUS	\$825
	<input type="checkbox"/> Prism Basic Back Bariatric Custom on chair -----	BBBCUS	\$1,075
	<input type="checkbox"/> Prism Ultra Back on chair Stock -----	UBOC	\$250
	<input type="checkbox"/> Prism Ultra Back on chair Non Stock -----	UBOC	\$300
	<input type="checkbox"/> Prism Ultra Back on chair Bariatric -----	UBOC	\$350
	<input type="checkbox"/> Prism Ultra Back Custom on chair-----	UBCUS	\$1,010
	<input type="checkbox"/> Prism Ultra Back Bariatric Custom on chair -----	UBBCUS	\$1,260
	<input type="checkbox"/> Prism Air Back on chair Stock -----	ABOC	\$400
	<input type="checkbox"/> Prism Air Back on chair Non Stock -----	ABOC	\$475
	<input type="checkbox"/> Prism Air Back on chair Bariatric -----	ABOC	\$600
	<input type="checkbox"/> Prism Air Back Custom on chair-----	ABCUS	\$1,075
	<input type="checkbox"/> Prism Air Back Bariatric Custom on chair -----	ABBCUS	\$1,325
	<input type="checkbox"/> Prism TrueFitt Back on chair Stock -----	TFBOC	\$650
	<input type="checkbox"/> Prism TrueFitt Back on chair Non Stock -----	TFBOC	\$700
	<input type="checkbox"/> Prism TrueFitt Back on chair Bariatric -----	TFBOC	\$950
	<input type="checkbox"/> Prism TrueFitt Back Custom on chair-----	TFBCUS	\$1,025
	<input type="checkbox"/> Prism TrueFitt Back Bariatric Custom on chair-----	TFBBCUS	\$1,350
	<input type="checkbox"/> Prism Max Ultra Back on chair Stock-----	MUBOC	\$300
	<input type="checkbox"/> Prism Max Ultra Back on chair Non Stock-----	MUBOC	\$350
	<input type="checkbox"/> Prism Max Ultra Back on chair Bariatric-----	MUBOC	\$400
	<input type="checkbox"/> Prism Max Ultra Back Custom on chair-----	MUBCUS	\$1,095
	<input type="checkbox"/> Prism Max Ultra Back Bariatric Custom on chair -----	MUBBCUS	\$1,345
	<input type="checkbox"/> Prism Max Air Back on chair Stock-----	MABOC	\$450
	<input type="checkbox"/> Prism Max Air Back on chair Non Stock-----	MABOC	\$525
	<input type="checkbox"/> Prism Max Air Back on chair Bariatric-----	MABOC	\$650
	<input type="checkbox"/> Prism Max Air Back Custom on chair-----	MABCUS	\$1,160
	<input type="checkbox"/> Prism Max Air Back Bariatric Custom on chair -----	MABBCUS	\$1,410
	<input type="checkbox"/> Prism Max TruFitt Back on chair Stock-----	MTFBOC	\$700
	<input type="checkbox"/> Prism Max TruFitt Back on chair Non Stock-----	MTFBOC	\$750
	<input type="checkbox"/> Prism Max TruFitt Back on chair Bariatric-----	MTFBOC	\$1,000
	<input type="checkbox"/> Prism Max TruFitt Back Custom on chair-----	MTFBCUS	\$1,160
	<input type="checkbox"/> Prism Max TruFitt Back Bariatric Custom on chair---	MTFBBCUS	\$1,435
9	<input type="checkbox"/> Lunar Back Low No Lateral on chair	LBLNOC	\$400
	<input type="checkbox"/> Lunar Back Low Basic Lateral on chair	LBLBOC	\$475
10	<input type="checkbox"/> Lunar Back Mid No Lateral on chair	LBMNOC	\$490
	<input type="checkbox"/> Lunar Back Mid Basic Lateral on chair	LBMBOC	\$550
	<input type="checkbox"/> Lunar Back Mid Ultra Lateral on chair	LBMUOC	\$610
11	<input type="checkbox"/> Lunar Back Full No Lateral on chair	LBFNOC	\$510
	<input type="checkbox"/> Lunar Back Full Basic Lateral on chair	LBFBOC	\$570
	<input type="checkbox"/> Lunar Back Full Ultra Lateral on chair	LBFUOC	\$630
	<input type="checkbox"/> Lunar Back Custom on chair	LBCUSON	\$950
	<input type="checkbox"/> Omit Back -----	OB	less net \$30
8	<input type="checkbox"/> Specify Max Back Color: _____		
	<input type="checkbox"/> Specify Back Size W _____ x H _____		

FOOTNOTES:

- 5. seat pan is adjustable from 17" deep to 20" deep
- 6. back cane is adjustable from 16" wide to 24" wide
- 7. back cane is adjustable vertically 3" in 1" increments
- 8. Max Backs come in one of the following colours:
Black, Navy Blue, Royal Blue, Dark Green, or Burgundy
- * 9. Sizes 14W, 16W, 18W, 20W x 10H or 12H
- * 10. Sizes 14W, 16W, 18W, 20W x 14H, 16H or 18H
- * 11. Sizes 16W, 18W, 20W, 22W x 20H, 22H or 24H
- * (refer to lunar backs order form for specific sizes)

SEAT OPTIONS

<input type="checkbox"/> Seat Pan -----	SSP	Standard
<input type="checkbox"/> Prism Ideal Cushion - Stock Sizes -----	ICOC	Standard
<input type="checkbox"/> Prism Ideal Cushion - Non-Stock Sizes -----	ICOC	No Charge
<input type="checkbox"/> Prism Ideal Cushion Custom on chair-----	ICCUS	No Charge
<input type="checkbox"/> Prism Elite Cushion - Stock Sizes -----	ECOC	\$250
<input type="checkbox"/> Prism Elite Cushion - Non-Stock Sizes -----	ECOC	\$275
<input type="checkbox"/> Prism Elite Cushion Custom on chair-----	ECCUS	\$425
<input type="checkbox"/> Prism Supreme Cushion - Stock Sizes -----	SCOC	\$375
<input type="checkbox"/> Prism Supreme Cushion -Non-Stock Sizes -----	SCOC	\$400
<input type="checkbox"/> Prism Supreme Cushion Custom on chair-----	SCCUS	\$685
<input type="checkbox"/> Prism Supreme II Plus Cushion Stock on chair-----	SCPOC	\$375
<input type="checkbox"/> Prism Supreme II Plus Cushion Non-Stock on chair-----	SCPOC	\$400
<input type="checkbox"/> Prism Supreme II Plus Cushion Custom on chair-----	SCPCUS	\$685
<input type="checkbox"/> Prism Supreme Air Cushion - Stock Sizes -----	SACOC	\$500
<input type="checkbox"/> Prism Supreme Air Cushion -Non-Stock Sizes -----	SACOC	\$575
<input type="checkbox"/> Prism Supreme Air Cushion Custom on chair-----	SACCUS	\$785
<input type="checkbox"/> Omit Cushion -----	OC	less net \$10

WARRANTY
<p>Future Mobility Healthcare Inc. warrants the Orion Explorer Magellan as follows:</p> <p>a. Base frame, seat frame and wheel hubs (not tires) 60 months limited warranty</p> <p>b. Controller, charger and cables 24 months limited warranty</p> <p>c. Motors and gearboxes 18 months limited warranty</p> <p>d. Upholstery, plastic and rubber components, bearings, front rigging, forks and all other parts 6 months limited warranty</p>

ARMPADS

<input type="checkbox"/> Oval Full Length, Right-----	OFLR	Standard
<input type="checkbox"/> Oval Full Length, Left-----	OFLR	Standard
<input type="checkbox"/> Wide Multiposition Full Length, Right-----	WMFR	\$64
<input type="checkbox"/> Wide Multiposition Full Length, Left-----	WMFL	\$64

BACK ANGLE

<input type="checkbox"/> Fixed (+3 to -20 degrees)-----	FIXED	Standard
<input type="checkbox"/> Manual Recline (+5 to -40 degrees)-----	MANUAL	\$365

FRONT RIGGING

<input type="checkbox"/> 60 Degree Swingaway Upgraded-----	SAWAY60	\$150
<input type="checkbox"/> 70 Degree Swingaway Upgraded-----	SAWAY	\$150
<input type="checkbox"/> 90 Degree Swingaway Upgraded-----	SAWAY90	\$175
<input type="checkbox"/> 70 Degree Pin Style-----	70PIN	No Charge
<input type="checkbox"/> 80 Degree Pin Style-----	80PIN	No Charge
<input type="checkbox"/> Elevating Legrest (pair) -----	ELR	\$361
<input type="checkbox"/> Omit Front Rigging -----	OFR	less net \$20
<input type="checkbox"/> One Piece Footboard -----	OPFB	Standard

FOOTPLATE TYPE

<input type="checkbox"/> Composite Adjustable Angle Flip-Up Footplates-----	AAFPC	\$113
<input type="checkbox"/> Aluminum Adjustable Angle Flip-Up Footplates-----	AAFP	\$113
<input type="checkbox"/> Composite Flip-Up Footplates-----	EMCF	No Charge
<input type="checkbox"/> Adjustable Angle Flip-Up Spring Return Footplates -----	FFAASR	\$150

FOOTREST /LEGREST ACCESSORIES

<input type="checkbox"/> Heel Loops -----	HL	Standard
<input type="checkbox"/> Omit Heel Loops -----	OHL	less net \$7
<input type="checkbox"/> Foam Guard on Elevating Legrest (pair) -----	FG	\$70
<input type="checkbox"/> Foam Pad One Piece Footboard -----	OPFD	\$50
<input type="checkbox"/> Foam Pad on Footplate (pair) -----	FPF	\$70

HEADREST HARDWARE OPTIONS

<input type="checkbox"/> Fixed Headrest Hardware -----	OFH1	No Charge
<input type="checkbox"/> Single Multi-Axial Headrest Hardware -----	OFH2	\$158
<input type="checkbox"/> Dual Multi-Axial Headrest hardware -----	OFH3	\$195
<input type="checkbox"/> Fixed Headrest Hardware Ass'y (Plus 4")-----	OFH4	\$20
<input type="checkbox"/> Single Multi-Axial Headrest Hardware (Plus 4") -----	OFH5	\$183
<input type="checkbox"/> Dual Multi-Axial Headrest hardware (Plus 4") -----	OFH6	\$220

HEADREST FOAM PAD OPTIONS

<input type="checkbox"/> Small Headrest Foam Pad -----	OHP1	No Charge
<input type="checkbox"/> Small Headrest with Wings Foam Pad -----	OHP2	\$90
<input type="checkbox"/> Large Headrest Foam Pad -----	OHP3	No Charge
<input type="checkbox"/> Large Headrest with Wings Foam Pad -----	OHP4	\$95
<input type="checkbox"/> Omit Headrest -----	OHT	No Charge

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> Chest Belt-----	CBELT	\$106
<input type="checkbox"/> Anchor Brackets -----	ABKT	\$60
<input type="checkbox"/> Omit Lap Belt-----	OLB	No Charge