

Galaxy Lite

Order Form/Price List

(Effective August 6, 2010)

Account Information:

Account #: _____	PO# _____
Name: _____	Ship To: _____
Address: _____	Address: _____
City: _____	City: _____
Province: _____	Province: _____
Postal: _____	Postal: _____
Phone: _____	Phone: _____
Fax: _____	Fax: _____
Comment: _____	

1 FRAME TYPE

Folding Frame w/Swingaway Front End GALAXY \$2,350

FRAME FINISH

Black ----- FBK No Charge
 Silver ----- FSR No Charge
 Custom Colours ----- FCUS NET \$150

SEAT WIDTH

14" Seat Width ----- SW14 No Charge
 15" Seat Width ----- SW15 No Charge
 16" Seat Width ----- SW16 No Charge
 17" Seat Width ----- SW17 No Charge
 18" Seat Width ----- SW18 No Charge
 19" Seat Width ----- SW19 No Charge
 20" Seat Width ----- SW20 No Charge

SEAT DEPTH

15" Seat Depth ----- SD15 No Charge
 16" Seat Depth ----- SD16 No Charge
 17" Seat Depth ----- SD17 No Charge
 18" Seat Depth ----- SD18 No Charge
 19" Seat Depth ----- SD19 No Charge
 20" Seat Depth ----- SD20 No Charge

2 SEAT-TO-FLOOR HEIGHT - LOW FRAME (HEMI)

3	<input type="checkbox"/> 13.0" (20" rear & 4" caster & no spacer) -----	SFG13	\$225
4	<input type="checkbox"/> 14.0" (20" rear & 4" caster & no spacer) -----	SFG01	No Charge
	<input type="checkbox"/> 14.5" (20" rear & 4" caster & no spacer) -----	SFG02	No Charge
	<input type="checkbox"/> 15.0" (20" rear & 4" caster & no spacer) -----	SFG03	No Charge
	<input type="checkbox"/> 15.0" (20" rear & 5" caster & no spacer) -----	SFG05	No Charge
	<input type="checkbox"/> 15.5" (20" rear & 5" caster & no spacer) -----	SFG06	No Charge
5	<input type="checkbox"/> 15.5" (22" rear & 4" caster & 1" spacer) -----	SFG14	No Charge
	<input type="checkbox"/> 16.0" (20" rear & 4" caster & 1" spacer) -----	SFG17	No Charge
	<input type="checkbox"/> 16.0" (22" rear & 5" caster & 1" spacer) -----	SFG18	No Charge
	<input type="checkbox"/> 16.5" (20" rear & 5" caster & 1" spacer) -----	SFG19	No Charge

2 SEAT-TO-FLOOR HEIGHT - STANDARD FRAME

5	<input type="checkbox"/> 17.0" (20" rear & 6" caster) -----	SFG20	No Charge
	<input type="checkbox"/> 17.0" (22" rear & 5" caster & no spacer) -----	SFG07	No Charge
	<input type="checkbox"/> 17.5" (20" rear & 6" caster) -----	SFG21	No Charge
	<input type="checkbox"/> 17.5" (22" rear & 5" caster & 1" spacer) -----	SFG08	No Charge
	<input type="checkbox"/> 18.0" (20" rear & 6" caster) -----	SFG22	No Charge
	<input type="checkbox"/> 18.0" (22" rear & 6" caster & 1" spacer) -----	SFG09	No Charge
	<input type="checkbox"/> 18.0" (24" rear & 6" caster) -----	SFG23	No Charge
	<input type="checkbox"/> 18.5" (22" rear & 6" caster & 1" spacer) -----	SFG10	No Charge
	<input type="checkbox"/> 18.5" (24" rear & 6" caster) -----	SFG24	No Charge
	<input type="checkbox"/> 19.0" (22" rear & 7" caster & 1" spacer) -----	SFG12	No Charge
	<input type="checkbox"/> 19.5" (22" rear & 7" caster & 1" spacer) -----	SFG15	No Charge

Basic Chair Includes:

Lightweight Folding Frame
 Elliptical Heavyduty Crossbrace
 Tubing
 Seat Width (14"-20")
 Seat Depth (15"-20")
 Seat-To-Floor Height Front
 (13" - 19-1/2")
 Seat-To-Floor Height Rear
 (12-1/2" -18-1/4")
 Rear Wheel Options
 (20", 22" or 24")
 Handrim Options (Stainless Steel
 or Plastic Coated)
 Caster Options (4",5",6" or 7")
 Slip on Seat Upholstery
 Arm rest height adjustable
 1/2" increments
 Slip On Back Upholstery
 Angle Adjustable Push Canes
 Aluminum Flip-Up Footplate
 Heel Loops
 Wheel Locks (Push or Pull to Lock)
 Choice of Two Colors (Black
 or Silver Powdercoat Finish)

FOOTNOTES:

1. Seat angle is set at 3 degrees.
See service manual for details
2. Can not use a 8" caster with this wheelchair
3. Short fork is used to achieve these front seat to floor heights
4. For front seat to floor of 14" and above the long fork is used
5. Spacers are used in the front in order to maintain the seat angle at 3°

WHEEL OPTIONS

<input type="checkbox"/> Airless -----	WHLA	Standard
<input type="checkbox"/> Pneumatic -----	WHLP	No Charge
<input type="checkbox"/> Pneumatic with spokes -----	WHLPS	No Charge

CAMBER OPTIONS

<input type="checkbox"/> 0° Camber -----	CAMN0	Standard
<input type="checkbox"/> 3° Camber -----	CAMN3	No Charge

* 6. Sizes 14W, 16W, 18W, 20W x 10H or 12H
 * 7. Sizes 14W, 16W, 18W, 20W x 14H, 16H or 18H
 * 8. Sizes 16W, 18W, 20W, 22W x 20H, 22H or 24H
 * (refer to lunar backs order form for specific sizes)

HANDRIM OPTIONS

<input type="checkbox"/> 20" Stainless Steel -----	HSG20	Standard
<input type="checkbox"/> 22" Stainless Steel -----	HSG22	Standard
<input type="checkbox"/> 24" Stainless Steel -----	HSG24	Standard
<input type="checkbox"/> 20" Plastic Coated -----	HPG20	\$110
<input type="checkbox"/> 22" Plastic Coated -----	HPG22	\$110
<input type="checkbox"/> 24" Plastic Coated -----	HPG24	\$110
<input type="checkbox"/> Omit Handrims -----	OHR	less net \$15

ARM TYPE

<input type="checkbox"/> Full Length Adj. Height (vinyl pad) T-Style -----	AFLG	\$215
<input type="checkbox"/> Desk Length Adj. Height (vinyl pad) T-Style -----	ADLG	\$210
<input type="checkbox"/> Curved Gel Adj. Height T-Style -----	ACGG	\$275
<input type="checkbox"/> Curved Foam Adj. Height T-Style -----	ACFG	\$265
<input type="checkbox"/> Cantilever Arm Full Length Adj. Angle (vinyl pad) -----	ACFLK	\$225
<input type="checkbox"/> Cantilever Arm Desk Length Adj. Angle (vinyl pad) -----	ACDLK	\$225

BACK TYPE

<input type="checkbox"/> Angle Adj 10 deg Bend w/Push Handle (low) -----	BLPG	Standard
<input type="checkbox"/> Angle Adj 10 deg Bend w/Push Handle (high) -----	BHPG	\$90

BACK UPHOLSTERY

<input type="checkbox"/> Slip-on Padded Back Upholstery -----	SOG	Standard
<input type="checkbox"/> Adjustable Tension Back -----	TBG	\$175
<input type="checkbox"/> Prism Basic Back on chair Stock -----	SBBOC	\$500
<input type="checkbox"/> Prism Basic Back on chair Non Stock -----	SBBOC	\$525
<input type="checkbox"/> Prism Basic Back on chair Custom -----	SBBCUS	\$825
<input type="checkbox"/> Prism Ultra Back on chair Stock -----	SUBOC	\$600
<input type="checkbox"/> Prism Ultra Back on chair Non Stock -----	SUBOC	\$625
<input type="checkbox"/> Prism Ultra Back on chair Custom -----	SUBCUS	\$1,010
<input type="checkbox"/> Prism Air Back on chair Stock -----	SABOC	\$650
<input type="checkbox"/> Prism Air Back on chair Non Stock -----	SABOC	\$675
<input type="checkbox"/> Prism Air Back on chair Custom -----	SABCUS	\$1,075
<input type="checkbox"/> Prism TrueFitt Back on chair Stock -----	STFBOC	\$750
<input type="checkbox"/> Prism TrueFitt Back on chair Non Stock -----	STFBOC	\$775
<input type="checkbox"/> Prism TrueFitt Back on chair Custom -----	STFBCUS	\$1,025
<input type="checkbox"/> Prism Max Ultra Back on chair Stock -----	SMUBOC	\$650
<input type="checkbox"/> Prism Max Ultra Back on chair Non Stock -----	SMUBOC	\$675
<input type="checkbox"/> Prism Max Ultra Back Custom on chair -----	SMUBCUS	\$1,060
<input type="checkbox"/> Prism Max Air Back on chair Stock -----	SMABOC	\$700
<input type="checkbox"/> Prism Max Air Back on chair Non Stock -----	SMABOC	\$725
<input type="checkbox"/> Prism Max Air Back Custom on chair -----	SMABCUS	\$1,125
<input type="checkbox"/> Prism Max TruFitt Back on chair Stock -----	SMTFBOC	\$800
<input type="checkbox"/> Prism Max TruFitt Back on chair Non Stock -----	SMTFBOC	\$825
<input type="checkbox"/> Prism Max TruFitt Back Custom on chair -----	SMTFBCUS	\$1,075
6 <input type="checkbox"/> Lunar Back Low No Lateral on chair -----	SLBLNOC	\$500
<input type="checkbox"/> Lunar Back Low Basic Lateral on chair -----	SLBLBOC	\$575
7 <input type="checkbox"/> Lunar Back Mid No Lateral on chair -----	SLBMNOC	\$590
<input type="checkbox"/> Lunar Back Mid Basic Lateral on chair -----	SLBMBOC	\$650
<input type="checkbox"/> Lunar Back Mid Ultra Lateral on chair -----	SLBMUOC	\$710
8 <input type="checkbox"/> Lunar Back Full No Lateral on chair -----	SLBFNOC	\$610
<input type="checkbox"/> Lunar Back Full Basic Lateral on chair -----	SLBFBOC	\$670
<input type="checkbox"/> Lunar Back Full Ultra Lateral on chair -----	SLBFUOC	\$730
<input type="checkbox"/> Lunar Back Custom on chair -----	SLBCUSON	\$1,050
<input type="checkbox"/> Omit Back Upholstery -----	OBU	less net \$20
<input type="checkbox"/> Specify Back Size W____x H____		

SEAT OPTIONS

<input type="checkbox"/> Black Nylon Seat Upholstery -----	SUG	Standard
<input type="checkbox"/> Prism Ideal Cushion - Stock Sizes -----	SICOC	\$175
<input type="checkbox"/> Prism Ideal Cushion - Non-Stock Sizes -----	SICOC	\$200
<input type="checkbox"/> Prism Ideal Cushion - Custom Sizes -----	ICCUS	\$275
<input type="checkbox"/> Prism Elite Cushion - Stock Sizes -----	ECOC	\$250
<input type="checkbox"/> Prism Elite Cushion - Non-Stock Sizes -----	ECOC	\$275
<input type="checkbox"/> Prism Elite Cushion - Custom Sizes -----	ECCUS	\$425
<input type="checkbox"/> Prism Supreme Cushion - Stock Sizes -----	SCOC	\$375
<input type="checkbox"/> Prism Supreme Cushion - Non-Stock Sizes -----	SCOC	\$400
<input type="checkbox"/> Prism Supreme Cushion - Custom Sizes -----	SCCUS	\$685
<input type="checkbox"/> Prism Supreme II Plus Cushion Stock -----	SCPOC	\$375
<input type="checkbox"/> Prism Supreme II Plus Cushion Non-Stock -----	SCPOC	\$400
<input type="checkbox"/> Prism Supreme II Plus Cushion Custom -----	SCPCUS	\$685
<input type="checkbox"/> Prism Supreme Air Cushion - Stock Sizes -----	SACOC	\$500
<input type="checkbox"/> Prism Supreme Air Cushion - Non-Stock Sizes -----	SACOC	\$575
<input type="checkbox"/> Prism Supreme Air Cushion - Custom Sizes -----	SACCUS	\$785
<input type="checkbox"/> Prism Supreme Drop Base Flushmount Stock Sizes----	SDB0OC	\$755
<input type="checkbox"/> Prism Supreme Drop Base 1" Stock Sizes -----	SDB1OC	\$755
<input type="checkbox"/> Prism Supreme Drop Base 2" Stock Sizes -----	SDB2OC	\$755
<input type="checkbox"/> Prism Supreme Drop Base Flushmount Non-Stock ----	SDB0OC	\$810
<input type="checkbox"/> Prism Supreme Drop Base 1" Non Stock -----	SDB1OC	\$810
<input type="checkbox"/> Prism Supreme Drop Base 2" Non Stock -----	SDB2OC	\$810
<input type="checkbox"/> Prism Supreme Drop Base Flushmount Custom -----	SDB0CUS	\$935
<input type="checkbox"/> Prism Supreme Drop Base 1" Custom Sizes -----	SDB1CUS	\$935
<input type="checkbox"/> Prism Supreme Drop Base 2" Custom Sizes -----	SDB2CUS	\$935
<input type="checkbox"/> Prism Drop Base Flushmount Stock Sizes -----	DB0OC	\$180
<input type="checkbox"/> Prism Drop Base 1" Stock Sizes -----	DB1OC	\$180
<input type="checkbox"/> Prism Drop Base 2" Stock Sizes -----	DB2OC	\$180
<input type="checkbox"/> Prism Drop Base Flushmount Non Stock Sizes-----	DB0OC	\$210
<input type="checkbox"/> Prism Drop Base 1" Non Stock Sizes -----	DB1OC	\$210
<input type="checkbox"/> Prism Drop Base 2" Non Stock Sizes -----	DB2OC	\$210
<input type="checkbox"/> Prism Drop Base Flushmount Custom Sizes -----	DB0CUS	\$250
<input type="checkbox"/> Prism Drop Base 1" Custom Sizes -----	DB1CUS	\$250
<input type="checkbox"/> Prism Drop Base 2" Custom Sizes -----	DB2CUS	\$250
<input type="checkbox"/> Omit Seat Upholstery -----	OSU	less net \$20

FRONT RIGGING

<input type="checkbox"/> 60 Degree Swingaway Upgraded -----	A60G	\$210
<input type="checkbox"/> 70 Degree Swingaway Upgraded -----	A70G	\$210

FOOTPLATE TYPE

<input type="checkbox"/> Composite Flip-Up Footplates -----	CFG	No Charge
<input type="checkbox"/> Aluminum Flip-Up Footplates -----	AFG	Standard
<input type="checkbox"/> X-Large Aluminum Flip-Up Footplates -----	XAFG	\$90
<input type="checkbox"/> Composite Adjustable Angle Flip-Up Footplates -----	FFAAC	\$120
<input type="checkbox"/> Aluminum Adjustable Angle Flip-Up Footplates -----	FFAA	\$120
<input type="checkbox"/> Composite Adjustable Angle Flip-Up Footplates Spring Return	FFAASR	\$150
<input type="checkbox"/> One Piece Flip-Up Footboard -----	OPFG	\$160

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> Padded Calf Strap -----	CSG	\$70
<input type="checkbox"/> Foam Pad on Footplate (pair) -----	FPF	\$70
<input type="checkbox"/> Foam Pad on One Piece Footboard -----	OPFD	\$50
<input type="checkbox"/> Heel Loops -----	HLG	Standard
<input type="checkbox"/> Omit Heel Loops -----	OHL	less net \$7

WHEEL LOCKS

<input type="checkbox"/> Push-to-Lock Wheel Lock -----	PSG	Standard
<input type="checkbox"/> Pull-to-Lock Wheel Lock -----	PLG	No Charge
<input type="checkbox"/> Wheel Lock Extension -----	WLG	No Charge

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> Anti-tippers -----	ANG	No Charge
<input type="checkbox"/> Auto-Seat Positioning Straps -----	ASG	\$52
<input type="checkbox"/> Velcro Release Positioning Strap -----	VRG	\$48
<input type="checkbox"/> Quick Release Axles -----	QRG	No Charge
<input type="checkbox"/> Clear Lap Tray 15"-18" -----	CLTSG	\$255
<input type="checkbox"/> Clear Lap Tray 19"-22" -----	CLTLG	\$255
<input type="checkbox"/> Padded Lap Tray 15"-16" -----	PLTG15	\$230
<input type="checkbox"/> Padded Lap Tray 17"-18" -----	PLTG17	\$230

- Padded Lap Tray 19"-20" ----- PLTG19 \$230
- Half Padded Lap Tray, left side ----- HPLTLG \$160
- Half Padded Lap Tray, right side ----- HPLTRG \$160

STF	REAR	CASTER	FORK	FORK	AXLE PLATE	SPACER	STF	STEM SIZE	SIDE
FRONT	WHEEL	SIZE	SIZE	HOLE	HOLE	SIZE	rear		FRAME
(inches)	(inches)	(inches)				(inches)	(inches)		STYLE
13.00	20	4	SHORT	A	SUPER LOW	0	12.50	SHORT	LOW
14.00	20	4	LRG	A	1	0	13.50	SHORT	LOW
14.50	20	4	LRG	B	1	0	14.00	SHORT	LOW
15.00	20	4	LRG	C	2	0	14.25	SHORT	LOW
15.00	20	5	LRG	B	2	0	14.25	SHORT	LOW
15.50	20	5	LRG	C	2	0	14.75	SHORT	LOW
15.50	22	4	LRG	B	2	1	14.75	LONG	LOW
16.00	20	4	LRG	C	3	1	15.25	LONG	LOW
16.00	22	5	LRG	B	2	1	15.25	LONG	LOW
16.50	20	5	LRG	C	3	1	16.00	LONG	LOW
17.00	20	6	LRG	B	3		16.25	SHORT	STANDARD
17.00	22	5	LRG	C	2	0	15.75	SHORT	STANDARD
17.50	20	6	LRG	C	4		17.00	SHORT	STANDARD
17.50	22	5	LRG	B	3	1	17.00	LONG	STANDARD
18.00	20	6	LRG	B	4		17.25	LONG	STANDARD
18.00	22	6	LRG	B	3	1	17.00	LONG	STANDARD
18.00	24	6	LRG	B	3		17.50	LONG	STANDARD
18.50	22	6	LRG	C	4	1	18.25	LONG	STANDARD
18.50	24	6	LRG	C	3		18.00	LONG	STANDARD
19.00	22	7	LRG	C	4	1	18.25	LONG	STANDARD
19.50	22	7	LRG	C	4	1+ ADJUST	18.25	LONG	STANDARD

