

## Headrest Options

**Fixed  
Single Multi-Axial  
Dual Multi-Axial**



### **Fixed Hardware**

- ◆ Horizontal and vertical adjustable hardware
- ◆ Available with Large, Small or with Wings Headrest and foam pad



### **Single Multi-Axial Hardware**

- ◆ Horizontal, vertical and single multi-axial adjustable hardware
- ◆ Available with Large, Small or with Wings Headrest and foam pad



### **Dual Multi-Axial Hardware**

- ◆ Horizontal, vertical and dual multi-axial adjustable hardware
- ◆ Available with Large, Small or with Wings Headrest and foam pad



### **Small Headrest Foam Pad**

- ◆ Contoured head support
- ◆ Constructed of layered foam on a coated metal plate, covered with soft stretch material
- ◆ Removable, washable cover
- ◆ Size: 12"Wx4"H



### **Large Headrest Foam Pad**

- ◆ Contoured head support
- ◆ Constructed of layered foam on a coated metal plate, covered with soft stretch material
- ◆ Removable, washable cover
- ◆ Size: 12"Wx9"H



### **Small Headrest Foam Pad with Wings**

- ◆ Anatomically shaped for more comfort
- ◆ Constructed of layered foam on a curved coated metal plate, covered with soft stretch material
- ◆ Removable, washable cover
- ◆ Size: 14"Wx4"H



### **Large Headrest Foam Pad with Wings**

- ◆ Anatomically shaped for more comfort
- ◆ Constructed of layered foam on a curved coated metal plate, covered with soft stretch material
- ◆ Removable, washable cover
- ◆ Size: 14"Wx9"H