

Lap Tray



The Lap Tray was designed for clients who require additional support in keeping items in reach.

Ideal for eating, writing, working, doing crafts and playing games.

The ultra rigid plastic tray has a smooth flat top with a raised spill edge all around, while the padded tray offers comfort and padding with non-slip vinyl.



The Lap Tray can be attached to the Orion Wheelchair or with additional hardware it can be adapted to other chairs in the market.



SPOTLIGHT ON

- Two styles—clear or padded
- Flexible sizes
- Attractive design

Clear Lap Tray Sizes

Medium 15"-18"
Large 19"-22"

Padded Lap Tray Sizes
16" - 24"

WARRANTY

One (1) years limited.

Future Mobility Healthcare