



Future Mobility Healthcare manufactures and markets products that increase the freedom and safety of mobility challenged individuals.

Renowned for quality and service, we provide high quality and innovative mobility products to the home healthcare market.

These Comfort Chairs are the latest stylish lift chairs on the market. Their look complements any home décor with a superior design that provides exceptional strength.

These chairs are enriched with several benefits, which provide style and durability for the user, all at a great value.

All chairs come with a lifetime warranty on the inner metal frame and 3 year warranty on the Okin actuators and motors.

Colour Selection

Faux Leather:

- Black
- Brown
- Rust



Microfiber:

- Blue
- Burgundy
- Green
- Brown
- Cream



Lift Chairs Comfort Series

With standard
heat & massage

Separate
hand control

Heat on lower back

Memory Foam

Foot massage



Future Mobility Healthcare Inc.
3223 Orlando Drive
Mississauga, ON L4V 1C5
Tel: 905-671-1661 or 1-888-737-4011
Fax: 905-671-3377 or 1-866-670-3377
www.future-mobility.com
information@future-mobility.com

Future Mobility Healthcare
innovative seating and mobility solutions

AG-001



- Weight Capacity: 375 lbs
- Overall Height: 43"
- Overall Width: 38"
- Back Height from seat: 25.5"
- Width between arms: 19"
- Seat Depth: 20.5"
- Seat height from floor: 20"
- Reclined length, overall length: 68"
- Faux leather: Black & Brown

AG-002



- Weight Capacity: 375 lbs
- Overall Height: 41"
- Overall Width: 35"
- Back Height from seat: 27"
- Width between arms: 21"
- Seat Depth: 20"
- Seat height from floor: 19"
- Reclined length, overall length: 64"
- Faux leather: Black, Brown
- Microfiber: Blue, Burgundy, Brown & Green

STANDARD FEATURES

- 3 position lift (quiet & smooth)
- Separate control unit for heat & massage
- Memory Foam seat
- Okin Motor with battery back-up

AG-003



- Weight Capacity: 375 lbs
- Overall Height: 41"
- Overall Width: 33.5"
- Back Height from seat: 27"
- Width between arms: 18"
- Seat Depth: 20.5"
- Seat height from floor: 20"
- Reclined length, overall length: 67"
- Microfiber: Blue, Burgundy, Brown, Green & Cream

AG-004



- Weight Capacity: 375 lbs
- Overall Height: 38"
- Overall Width: 38"
- Back Height from seat: 26"
- Width between arms: 19"
- Seat Depth: 20.5"
- Seat height from floor: 20"
- Reclined length, overall length: 68"
- Microfiber: Blue, Burgundy, Brown & Green

AG-005



- Weight Capacity: 375 lbs
- Overall Height: 44"
- Overall Width: 35.5"
- Back Height from seat: 27"
- Width between arms: 19"
- Seat Depth: 20.5"
- Seat height from floor: 17"
- Reclined length, overall length: 67"
- Faux leather: Black & Brown

AG-006



- Weight Capacity: 375 lbs
- Overall Height: 38"
- Overall Width: 38"
- Back Height from seat: 26"
- Width between arms: 20"
- Seat Depth: 21"
- Seat height from floor: 20"
- Reclined length, overall length: 68"
- Faux leather: Black & Rust

See colour selection on the back of this brochure!