

New

MAX Back

Future Mobility Healthcare manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Renown for quality and service we bring innovation to all of our products for the home healthcare market.

The Prism™ MAX Backs were designed for users who require moderate to moderately complex postural support and offers the most aggressive lateral support in the Prism line of products.

Available as an upgrade to the Prism Ultra, Prism Air, and Prism TruFitt Backs, this type of back comes equipped with high density foam, 3" higher set laterals and a new cover made with a breathable air-mesh Spandex material, in a variety of 5 colours.

Combine this with Prism's sleek and attractive tailored look for a back that will appeal to a wide variety of users!





SPOTLIGHT ON

- Lightweight, sturdy aluminum shell
- Seat angle adjustment from 82 degrees to 114 degrees optimizing users' positioning.
- 3" higher set laterals than regular Prism backs
- Attractive cover design in a selection of black, royal blue, navy blue, dark green and burgundy.
- Cane clamps included
- Quick push button release for easy installation

WARRANTY

Three (3) years limited on Back and 90 days on cover, 1 year on hardware.

PRODUCT WEIGHT

5.3 lbs (7.3 lbs bariatric)


USER WEIGHT

250 lbs (350 lbs bariatric)



ADP code SEBFM0150

(an upcharge of \$100 in addition to the ADP codes for the Ultra, Air, or TruFitt Backs)



Stock sizes—same day shipments

16" x 16"	M_B1616	16" x 20"	M_B1620
18" x 16"	M_B1816	18" x 20"	M_B1820
20" x 16"	M_B2016	20" x 20"	M_B2020

Non-stock sizes—shipments within 3 business days

14" x 14"	M_B1414	18" x 18"	M_B1818
17" x 16"	M_B1716	17" x 20"	M_B1720
19" x 16"	M_B1916	19" x 20"	M_B1920
16" x 18"	M_B1618	20" x 18"	M_B2018

Custom sizes—shipments within 7 business days

M_BCUS

Bariatric sizes—shipments within 3 business days

21" x 16"	M_B2116	23" x 16"	M_B2316
21" x 18"	M_B2118	23" x 18"	M_B2318
21" x 20"	M_B2120	23" x 20"	M_B2320
22" x 16"	M_B2216	24" x 16"	M_B2416
22" x 18"	M_B2218	24" x 18"	M_B2418
22" x 20"	M_B2220	24" x 20"	M_B2420

Custom bariatric sizes—shipments within 7 business days

M_BBCUS

_ U for Max Ultra / A for Max Air / TF for Max TruFitt