

Northern Lite

Order Form/Price List

(Effective August 20, 2010)

Account Information:

Account #: _____
 Name: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____
 Comment: _____

Date: _____

PO# _____
 Ship To: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____

1 FRAME TYPE

Folding Frame w/Swingaway Front End NORTHLT \$1,650

FRAME FINISH

Silver FSR Standard

SEAT SIZES

16" Seat Width x 16" Seat Depth SW16D16 No Charge
 16" Seat Width x 18" Seat Depth SW16D18 No Charge
 18" Seat Width x 16" Seat Depth SW18D16 No Charge
 18" Seat Width x 18" Seat Depth SW18D18 No Charge
 20" Seat Width x 18" Seat Depth SW20D18 No Charge

FRONT SEAT-TO-FLOOR HEIGHT & CASTER SIZE - Choose one

13.75" (4" Caster) STF11 \$55
 14.5" (4" Caster) STF12 \$55
 15.75" (4" Caster) STF13 No Charge
 16.5" (4" Caster) STF14 No Charge
 16" (6" Caster) STF15 No Charge
 16.75" (6" Caster) STF16 Standard
 17.5" (6" Caster) STF17 No Charge
 17.75" (8" Caster) STF18 No Charge
 18.5" (8" Caster) STF19 No Charge
 19" (8" Caster) STF110 \$35

REAR SEAT-TO-FLOOR HEIGHT & WHEEL SIZE - Choose one

2 13" (20" Wheels) STF111 No Charge
 13.75" (20" Wheels) STF112 No Charge
 2 14" (22" Wheels) STF113 No Charge
 14.5" (20" Wheels) STF114 No Charge
 14.75" (22" Wheels) STF115 No Charge
 15.25" (20" Wheels) STF116 No Charge
 15.5" (22" Wheels) STF117 No Charge
 15.75" (24" Wheels) STF118 No Charge
 16.25" (22" Wheels) STF119 No Charge
 16.5" (24" Wheels) STF120 Standard
 17.25" (24" Wheels) STF121 No Charge
 18" (24" Wheels) STF122 No Charge
 18.75 (24" Wheels) STF123 No Charge

WHEEL OPTIONS

3 Airless (20" , 22" , 24") WHLNA Standard
 Pneumatic (24") WHLNP No Charge

HANDRIM OPTIONS

Aluminum HRA Standard
 Omit Handrims OHR less net \$15

Basic Chair Includes:

Lightweight Folding Frame w/Swingaway Front End
 Elliptical Heavyduty Double Crossbrace Tubing
 Seat Width (16"-20")
 Seat Depth (16"-18")
 Front Seat-To-Floor Height (14"-19")
 Rear Wheel Options (20", 22" or 24") Airless
 Handrim Options (Aluminum)
 Caster Options (4", 6" or 8")
 Camber Adjustable Rear Wheels (0°, 3°, 5°)
 Seat Upholstery
 Northern Lite Seat Cushion
 Adjustable Tension Back Upholstery
 Push Canes
 Composite Flip-Up Footplates
 Heel Loops
 Wheel Locks (Push or Pull to Lock)
 Full Length Adj. Height T-Arm (Vinyl Pad)
 Anti tippers
 Seat Upholstery Bottom Storage Pocket
 Silver Powdercoat Finish

FOOTNOTES:

1. Front seat height should be larger than rear seat height
2. Can only have 0° camber angle on rear wheels at these heights
3. 20" and 22" rear wheels available in airless only. 24" rear wheels available in both airless and pneumatic

ARM TYPE

- Full Length Adj. Height (Vinyl Pad) ----- ARMNFL Standard
- Desk Length Adj. Height (Vinyl Pad) ----- ARMNDL No Charge

BACK TYPE

- 4 Back 10 deg bend w/Push Handle (low) ----- BACKNL Standard
- 5 Back 10 deg bend w/Push Handle (high) ----- BACKNH No Charge

BACK UPHOLSTERY

- Adjustable Tension Back Upholstery ----- BUN Standard
- Prism Basic Back on chair Stock ----- SBBOC \$500
- Prism Basic Back on chair Non Stock ----- SBBOC \$525
- Prism Basic Back on chair Cus ----- SBBCUS \$825
- Prism Ultra Back on chair Stock ----- SUBOC \$600
- Prism Ultra Back on chair Non Stock ----- SUBOC \$625
- Prism Ultra Back on chair Cus ----- SUBCUS \$1,010
- Prism Air Back on chair Stock ----- SABOC \$650
- Prism Air Back on chair Non Stock ----- SABOC \$675
- Prism Air Back on chair Cus ----- SABCUS \$1,075
- Prism TrueFitt Back on chair Stock ----- STFBOC \$750
- Prism TrueFitt Back on chair Non Stock ----- STFBOC \$775
- Prism TrueFitt Back on chair Cus ----- STFBCUS \$1,025
- Prism Max Ultra Back on chair Stock ----- SMUBOC \$650
- Prism Max Ultra Back on chair Non Stock ----- SMUBOC \$675
- Prism Max Ultra Back on chair Cus ----- SMUBCUS \$1,060
- Prism Max Air Back on chair Stock ----- SMABOC \$700
- Prism Max Air Back on chair Non Stock ----- SMABOC \$725
- Prism Max Air Back on chair Cus ----- SMABCUS \$1,125
- Prism Max TrueFitt Back on chair Stock ----- SMTFBOC \$800
- Prism Max TrueFitt Back on chair Non Stock ----- SMTFBOC \$825
- Prism Max TrueFitt Back on chair Cus ----- SMTFBCUS \$1,075
- Omit Back Upholstery ----- OBU less net \$20
- Specify Back Size W _____ x H _____

SEAT OPTIONS

- Black Nylon Seat Upholstery ----- SUNL Standard
- Northern Lite Seat Cushion ----- NLC Standard
- Omit Seat Upholstery ----- OSU less net \$20
- Prism Ideal Cushion - Stock Sizes ----- SICOC \$175
- Prism Ideal Cushion - Non-Stock Sizes ----- SICOC \$200
- Prism Ideal Cushion - Custom Sizes ----- ICCUS \$275
- Prism Elite Cushion - Stock Sizes ----- ECOC \$250
- Prism Elite Cushion - Non-Stock Sizes ----- ECOC \$275
- Prism Elite Cushion - Custom Sizes ----- ECCUS \$425
- Prism Supreme Cushion - Stock Sizes ----- SCOC \$375
- Prism Supreme Cushion - Non-Stock Sizes ----- SCOC \$400
- Prism Supreme Cushion - Custom Sizes ----- SCCUS \$685
- Prism Supreme II Plus Cushion Stock ----- SCPOC \$375
- Prism Supreme II Plus Cushion Non-Stock ----- SCPOC \$400
- Prism Supreme II Plus Cushion Custom ----- SCPCUS \$685
- Prism Supreme Air Cushion - Stock Sizes ----- SACOC \$500
- Prism Supreme Air Cushion - Non-Stock Sizes ----- SACOC \$575
- Prism Supreme Drop Base Flushmount Stock ----- SDB0OC \$755
- Prism Supreme Drop Base 1" Stock Sizes ----- SDB1OC \$755
- Prism Supreme Drop Base 2" Stock Sizes ----- SDB2OC \$755
- Prism Supreme Drop Base Flushmount Non-Stock ----- SDB0OC \$810
- Prism Supreme Drop Base 1" Non Stock ----- SDB1OC \$810
- Prism Supreme Drop Base 2" Non Stock ----- SDB2OC \$810
- Prism Supreme Drop Base Flushmount Custom ----- SDB0CUS \$935
- Prism Supreme Drop Base 1" Custom ----- SDB1CUS \$935
- Prism Supreme Drop Base 2" Custom ----- SDB2CUS \$935
- Prism Drop Base Flushmount Stock Sizes ----- DB0OC \$180
- Prism Drop Base 1" Stock Sizes ----- DB1OC \$180
- Prism Drop Base 2" Stock Sizes ----- DB2OC \$180
- Prism Drop Base Flushmount Non-Stock ----- DB0OC \$210
- Prism Drop Base 1" Non Stock Sizes ----- DB1NS \$210
- Prism Drop Base 2" Non Stock Sizes ----- DB2NS \$210
- Prism Drop Base Flushmount Custom ----- DB0CUS \$250
- Prism Drop Base 1" Custom ----- DB1CUS \$250
- Prism Drop Base 2" Custom ----- DB2CUS \$250

FOOTNOTES:

- 4. Top of back cane is set at 16" high from the top of the seat rail
- 5. Top of back cane is set at 20" high from the top of the seat rail

FRONT RIGGING

- 70 Degree Swingaway Pinstyle ----- SW70N Standard
- Elevating Legrest (pair) ----- ELRN \$285

FOOTPLATE TYPE

- Composite Flip-Up Footplates ----- FPCN Standard
- Adjustable Angle Flip-Up Footplates ----- FFAAN \$179

FOOTREST/LEGREST ACCESSORIES

- Padded Calf Strap ----- PCS \$70
- Foam Guards on Elevating Legrest (pair) ----- FG \$70
- Heel Loops ----- HL Standard
- Omit Heel Loops ----- OHL less net \$7

WHEEL LOCKS

- 6 Push-to-Lock Wheel Lock Standard Mount ----- WLPUSN Standard
- 7 Push-to-Lock Wheel Lock Transfer Mount ----- WLPUTN No Charge
- Wheel Lock Extension ----- WLEXTN \$55
- Pull-to-Lock ----- WLPLN No Charge

REAR WHEEL CAMBER ANGLE

- 8 No Camber ----- CAMN0 Standard
- 3 ° camber ----- CAMN3 No Charge
- 5 ° camber ----- CAMN5 No Charge

MISCELLANEOUS ACCESSORIES

- Anti-tippers ----- ATS No Charge
- Auto-Seat Positioning Straps ----- ARB \$52
- Velcro Release Positioning Strap ----- VRB \$48

FOOTNOTES:

- 6. Standard Mount - Footrests swing further back to rear wheels
- 7. Transfer Mount - Brakes levers are below seat level when locked
- 8. Standard Wheelchairs are set at a 0 ° camber angle on rear wheels

Future Mobility Healthcare Inc. reserves the right to change prices and specifications without notice.

Note: All specifications and dimensions are approximate.
Prism and Orion Products are trademarks of Future Mobility Healthcare Inc.