

## Manual Wheelchairs



*Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market. Future Mobility Healthcare Inc. is a Canadian Manufacturer of Class I Medical Devices.*



Orion // is the most innovative and affordable tilt-in-space wheelchair on the market today. Developed for those individuals requiring adjustability to redistribute their body weight for optimal positioning purposes, the Orion // assists in allowing the user to function independently and to become more involved in their lifestyles.

With a true 45 degrees of tilt, optional 30 degrees of recline, standard on-chair adjustable seating system, incredible range of seat-to-floor heights and creative streamlined design.

The Orion //...dynamic, progressive, sure to become the most popular tilt chair under the stars!

**Future M  obility Healthcare**



### SPOTLIGHT ON:

- Push bar AND stroller handles
- Standard seating system consisting of: two-way contoured, padded headrest, Orion Back with 1" depth adjustability and Prism Ideal Cushion
- Optional locking collar stops on tilt and recline
- "anti-rattle" anti-tippers
- Quick adjustable contoured armrest pads
- Quick removable arm
- Arm location adjustable by 1 1/2" to allow for tray placement
- Multi-adjustable axle plate
- Spring return push-to-lock wheel lock
- Incredible low 13 1/4" seat-to-floor height
- Unique push button release on optional front rigging
- Innovative elliptical tubing
- Clean, sleek appearance



### FEATURES:

- 45 degrees of infinite tilt
- 30 degrees of infinite recline
- Seat widths: 15"-20" (350lbs)
- **Bariatric widths: 26" - 32" (700lbs)**
- Seat depths: 15"-20"
- Seat-to-floor heights: 13 1/4" - 19"
- Arm heights: 9 1/2" - 14" (adjustable)
- Front rigging: 70 degree standard
- Adjustable fold down back
- Fixed Back
- 12", 20", 22" or 24" rear wheels
- 4", 5", 7" or 8" casters
- Low or tall push bars and back heights
- Aluminum seat pan

### SPECIFICATIONS

Tilt Range: 3-48 degrees  
Recline Range: 6-36 degrees

### OVERALL DIMENSIONS

(18" x 16" tilt & recline, 16 1/2" STF)  
Chair width: 27 1/2"  
Chair length: 34 1/2"  
Chair height: 44"  
Chair weight: 65 lbs. (including seating system)

### WARRANTY

Frame is warranted for the lifetime of the original purchaser or user.

**Future Mobility Healthcare**

3223 Orlando Drive, Mississauga, ON, L4V 1C5, 905-671-1661, 888-737-4011, fax 866-670-3377 [www.future-mobility.com](http://www.future-mobility.com)  
Future Mobility Healthcare Inc. Orion and Prism are trademarks of Future Mobility Healthcare Inc.