



Orion II Order Form/Price List (Effective August 6, 2010)

Account Information:

Date: _____

Account #: _____
 Name: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____
 Comment: _____

PO# _____
 Ship To: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____

FRAME TYPE

Orion II ----- ORION **\$2,700**

FRAME FINISH

Black ----- FBK No Charge
 Silver ----- FSR No Charge
 Custom Colours ----- FCUS NET \$150

SEAT WIDTH

15" Seat Width ----- SW15 No Charge
 16" Seat Width ----- SW16 No Charge
 17" Seat Width ----- SW17 No Charge
 18" Seat Width ----- SW18 No Charge
 19" Seat Width ----- SW19 No Charge
 20" Seat Width ----- SW20 No Charge
¹ 21" Seat Width (Weight Limit 350 lbs) ----- SW21 \$475
¹ 22" Seat Width (Weight Limit 350 lbs) ----- SW22 \$475
^{1,2} 23" Seat Width (Weight Limit 350 lbs) ----- SW23 \$685
^{1,2} 24" Seat Width (Weight Limit 350 lbs) ----- SW24 \$685

SEAT DEPTH

15" Seat Depth ----- SD15 No Charge
 16" Seat Depth ----- SD16 No Charge
 17" Seat Depth ----- SD17 No Charge
 18" Seat Depth ----- SD18 No Charge
 19" Seat Depth ----- SD19 No Charge
 20" Seat Depth ----- SD20 \$110
 21" Seat Depth ----- SD21 \$225
 22" Seat Depth ----- SD22 \$225
 23" Seat Depth ----- SD23 \$225
 24" Seat Depth ----- SD24 \$225

SEAT-TO-FLOOR HEIGHT (up to 350 lbs weight capacity) - Orion II

⁴ 13 1/4" (20" rear & 4" caster) ----- STF204U \$225
 14 1/2" (20" rear & 5" caster) ----- STF205U No Charge
 15" (20" rear & 5" caster) ----- STF205M No Charge
 15 1/2" (22" rear & 5" caster) ----- STF225U No Charge
⁵ 15 1/2" (12" rear & 6" caster) ----- STF126S \$195
 16" (20" rear & 5" caster) ----- STF205L No Charge
⁵ 16" (12" rear & 7" caster) ----- STF127S \$195
 16" (24" rear & 6" caster) ----- STF246U No Charge
 16" (22" rear & 6" caster) ----- STF226M No Charge
 16 1/2" (24" rear & 6" caster) ----- STF246M No Charge
 16 1/2" (24" rear & 7" caster) ----- STF247U No Charge
 16 1/2" (22" rear & 7" caster) ----- STF227M No Charge
 17" (24" rear & 7" caster) ----- STF247M No Charge
 17" (22" rear & 8" caster) ----- STF228L No Charge
 18" (24" rear & 8" caster) ----- STF248L Standard
 19" (24" rear & 8" caster, 1" spacer) ----- STF248L1 No Charge

Standard Features:
Tilt Function (45 Degrees)
Seat Width 15"-22"
Seat Depth 15"-19"
Canes and Stroller Bar
Adjustable Height Full Length Arm
T-Style
Height of 9 1/2" - 14"
Adjustable Arm location for
Tray placement
Low Fixed Back
Seat to Floor 18"
Rear Wheels 24"
Casters 8"
Multi Adjustable Axel Plate
70 Degree Front Rigging w/footplate
Heel Loops
Spring Back Push to Lock Wheel Locks
Aluminum Seat Pan
Prism Ideal Cushion
Adjustable Orion Back
Two way Padded, Contoured Headrest
Anti-tippers
Choice of 2 colours - Black or Silver

FOOTNOTES:

1. Heavy Duty Package consisting of gusseted frame
2. Quick release axles are used on 23"-24" wide chairs.
4. Only available with pull to lock brakes.
5. With 12" rear wheel you can only push to lock brakes.

20" (24" rear & 8" caster, 2" spacer) ----- STF248L2 \$110

WHEEL OPTIONS

Airless ----- WHLA Standard
 Pneumatic ----- WHLP No Charge

HANDRIM OPTIONS

Stainless Steel ----- SS Standard
 Plastic Coated ----- PC \$160
 Omit Handrims ----- OHR less net \$15

ARM TYPE

T-Style Adjustable Height - Gel ----- ADG Standard
 T-Style Adjustable Height - Foam ----- ADF No Charge
 T-Style Adjustable Height - Vinyl 14"----- AVTS \$125
 Flip Back Adjustable Height - Gel ----- ADFG \$125
 Flip Back Adjustable Height - Foam ----- ADFF \$125
 Flip Back Adjustable Height - Vinyl 14"----- AVFB \$125

PUSH OPTIONS

Canes & Stroller Bars ----- CPB Standard
 Strollers Bars ----- PB No Charge

BACKREST OPTIONS

Low Fixed ----- FBL Standard
 Tall Fixed ----- FBT No Charge
 Straight Push Cane Low Fixed ----- SPCL \$150
 Straight Push Cane High Fixed ----- SPCH \$150
 Special Height Push Cane/Stroller (___Ht) ----- SHPC \$180
 Low Adjustable Angle Fold Down ----- ABL \$250
 Tall Adjustable Angle Fold Down ----- ABT \$250
 Low Recline ----- RBL \$500
 Tall Recline----- RBT \$500

FOOTNOTES:

6. Height adjustability of 4" and depth adjustability of 1"

BACK OPTIONS

⁶ Orion Back (16", 18", 20") available height only----- OBK Standard
 Prism Basic Back on chair Stock ----- BBOC \$150
 Prism Basic Back on chair Non Stock ----- BBOC \$175
 Prism Basic Back on chair Bariatric ----- BBOC \$250
 Prism Basic Back Custom on chair----- BBCUS \$825
 Prism Basic Back Bariatric Custom on chair----- BBBCUS \$1,075
 Prism Ultra Back on chair Stock ----- UBOC \$250
 Prism Ultra Back on chair Non Stock ----- UBOC \$300
 Prism Ultra Back on chair Bariatric ----- UBOC \$350
 Prism Ultra Back Custom on chair----- UBCUS \$1,010
 Prism Ultra Back Bariatric Custom on chair----- UBBCUS \$1,260
 Prism Air Back on chair Stock ----- ABOC \$400
 Prism Air Back on chair Non Stock ----- ABOC \$475
 Prism Air Back on chair Bariatric ----- ABOC \$600
 Prism Air Back Custom on chair----- ABCUS \$1,075
 Prism Air Back Bariatric Custom on chair ----- ABBCUS \$1,325
 Prism TrueFitt Back on chair Stock ----- TFBOC \$650
 Prism TrueFitt Back on chair Non Stock ----- TFBOC \$700
 Prism TrueFitt Back on chair Bariatric ----- TFBOC \$950
 Prism TrueFitt Back Custom on chair----- TFBCUS \$1,025
 Prism TrueFitt Back Bariatric Custom on chair----- TFBBCUS \$1,350
 Prism Max Ultra Back on chair Stock----- MUBOC \$300
 Prism Max Ultra Back on chair Non Stock----- MUBOC \$350
 Prism Max Ultra Back on chair Bariatric----- MUBOC \$400
 Prism Max Ultra Back Custom on chair----- MUBCUS \$1,095
 Prism Max Ultra Back Bariatric Custom on chair ----- MUBBCUS \$1,345
 Prism Max Air Back on chair Stock----- MABOC \$450
 Prism Max Air Back on chair Non Stock----- MABOC \$525
 Prism Max Air Back on chair Bariatric----- MABOC \$650
 Prism Max Air Back Custom on chair----- MABCUS \$1,160

<input type="checkbox"/>	Prism Max Air Back Bariatric Custom on chair -----	MABBCUS	\$1,410
<input type="checkbox"/>	Prism Max TruFitt Back on chair Stock-----	MTFBOC	\$700
<input type="checkbox"/>	Prism Max TruFitt Back on chair Non Stock-----	MTFBOC	\$750
<input type="checkbox"/>	Prism Max TruFitt Back on chair Bariatric-----	MTFBOC	\$1,000
<input type="checkbox"/>	Prism Max TruFitt Back Custom on chair-----	MTFBCUS	\$1,160
<input type="checkbox"/>	Prism Max TruFitt Back Bariatric Custom on chair-----	MTFBBCUS	\$1,435
16	<input type="checkbox"/> Lunar Back Low No Lateral on chair -----	LBLNOC	\$400
<input type="checkbox"/>	Lunar Back Low Basic Lateral on chair -----	LBLBOC	\$475
17	<input type="checkbox"/> Lunar Back Mid No Lateral on chair -----	LBMNOC	\$490
<input type="checkbox"/>	Lunar Back Mid Basic Lateral on chair -----	LBMBOC	\$550
<input type="checkbox"/>	Lunar Back Mid Ultra Lateral on chair -----	LBMUOC	\$610
18	<input type="checkbox"/> Lunar Back Full No Lateral on chair -----	LBFNOC	\$510
<input type="checkbox"/>	Lunar Back Full Basic Lateral on chair -----	LBFBOC	\$570
<input type="checkbox"/>	Lunar Back Full Ultra Lateral on chair -----	LBFUOC	\$630
<input type="checkbox"/>	Lunar Back Custom on chair -----	LBCUSON	\$950
<input type="checkbox"/>	Omit Back -----	OB	less net \$30
7	<input type="checkbox"/> Specify Max Back Color: _____		
8	<input type="checkbox"/> Specify Back Size W ____ x H ____		

CUSTOM MODIFICATIONS

<input type="checkbox"/>	Longer Pins for Back Hardware 1/2"-----	BLPH	\$50
<input type="checkbox"/>	Longer Pins for Back Hardware 1"-----	BLPO	\$50
<input type="checkbox"/>	Longer Pins for Back Hardware 2"-----	BLPT	\$50

SEAT OPTIONS

<input type="checkbox"/>	Seat Pan -----	SSP	Standard
<input type="checkbox"/>	Prism Ideal Cushion - Stock Sizes -----	ICOC	Standard
<input type="checkbox"/>	Prism Ideal Cushion - Non-Stock Sizes -----	ICOC	No Charge
<input type="checkbox"/>	Prism Ideal Cushion Custom on chair-----	ICCUS	No Charge
<input type="checkbox"/>	Prism Elite Cushion - Stock Sizes -----	ECOC	\$250
<input type="checkbox"/>	Prism Elite Cushion - Non-Stock Sizes -----	ECOC	\$275
<input type="checkbox"/>	Prism Elite Cushion Custom on chair-----	ECCUS	\$425
<input type="checkbox"/>	Prism Supreme Cushion - Stock Sizes -----	SCOC	\$375
<input type="checkbox"/>	Prism Supreme Cushion -Non-Stock Sizes -----	SCOC	\$400
<input type="checkbox"/>	Prism Supreme Cushion Custom on chair-----	SCCUS	\$685
<input type="checkbox"/>	Prism Supreme II Plus Cushion - Stock Sizes -----	SCPOC	\$375
<input type="checkbox"/>	Prism Supreme II Plus Cushion - Non-Stock Sizes -----	SCPOC	\$400
<input type="checkbox"/>	Prism Supreme II Plus Cushion - Custom on Chair -----	SCPCUS	\$685
<input type="checkbox"/>	Prism Supreme Air Cushion - Stock Sizes -----	SACOC	\$500
<input type="checkbox"/>	Prism Supreme Air Cushion -Non-Stock Sizes -----	SACOC	\$575
<input type="checkbox"/>	Prism Supreme Air Cushion Custom on chair-----	SACCUS	\$785
<input type="checkbox"/>	Omit Cushion -----	OC	less net \$10

FRONT RIGGING

<input type="checkbox"/>	70 Degree Pin Style -----	FA70	Standard
<input type="checkbox"/>	80 Degree Pin Syle -----	FA80	\$135
<input type="checkbox"/>	60 Degree Swingaway Upgraded -----	F60	\$215
<input type="checkbox"/>	70 Degree Swingaway Upgraded -----	F70	\$215
<input type="checkbox"/>	90 Degree Swingaway Upgraded -----	F90	\$230
9	<input type="checkbox"/> Elevating Legrest (pair) -----	ELR	\$550
<input type="checkbox"/>	Omit Front Rigging -----	OFR	less net \$20
<input type="checkbox"/>	Omit Footrest Hanger -----	OFH	less net \$10

FOOTPLATE TYPE

<input type="checkbox"/>	Composite Flip Up Footplates-----	FPC	No Charge
<input type="checkbox"/>	Aluminum Flip Up Footplates-----	FPA	Standard
<input type="checkbox"/>	X-Large Aluminum Flip-Up Footplates -----	XFPA	\$90
<input type="checkbox"/>	X-Large Composite Flip-Up Footplates -----	XFPC	\$110
<input type="checkbox"/>	Composite Adjustable Angle Flip-Up Footplates -----	FFAAC	\$120
<input type="checkbox"/>	Aluminum Adjustable Angle Flip Up Footplates-----	FFAA	\$120
<input type="checkbox"/>	Composite Adjustable Angle Flip-Up Footplates Spring Return -----	FFAASR	\$150
10	<input type="checkbox"/> One Piece Footboard (15"-20")-----	FPOP	\$160
10	<input type="checkbox"/> One Piece Footboard (21"-24" Custom)-----	FPOPC	\$195

FOOTREST /LEGREST ACCESSORIES

<input type="checkbox"/>	Heel Loops -----	HL	Standard
<input type="checkbox"/>	Omit Heel Loops -----	OHL	less net \$7

FOOTNOTES:

- 7. Max Backs come in one of the following colours: Black, Navy Blue, Royal Blue, Forest Green or Burgundy
- 8. Please refer to Backs order form for for stock, non-stock, bariatric and custom sizes
- 9. Includes Swingway Calf Pads
- 10. Not available with 7" or 8" casters when footboard is facing back on 90 degree swingaway
- 11. Swingaway Calf Pad can only be assembled onto Swingaway front rigging.
- * 16. Sizes 14W, 16W, 18W, 20W x 10H or 12H
- * 17. Sizes 14W, 16W, 18W, 20W x 14H, 16H or 18H
- * 18. Sizes 16W, 18W, 20W, 22W x 20H, 22H or 24H
- * (refer to lunar backs order form for specific sizes)

<input type="checkbox"/>	Padded Calf Strap -----	PCS	\$70
11	<input type="checkbox"/>	Swingway Calf Pad (pair) -----	SCPP \$185
	<input type="checkbox"/>	Foam Guard on Elevating Legrest (pair) -----	FG \$70
	<input type="checkbox"/>	Foam Pad One Piece Footboard -----	OPFD \$50
	<input type="checkbox"/>	Foam Pad on Footplate (pair) -----	FPF \$70

WHEEL LOCKS

<input type="checkbox"/>	Push-to-Lock Wheel Lock -----	WLP	Standard
<input type="checkbox"/>	Pull-to-Lock Wheel Lock -----	WLPL	\$40
<input type="checkbox"/>	Wheel Lock Extension -----	WLE	No Charge

HEADREST HARDWARE OPTIONS

<input type="checkbox"/>	Fixed Headrest Hardware -----	OFH1	Standard
<input type="checkbox"/>	Single Multi-Axial Headrest Hardware -----	OFH2	\$158
<input type="checkbox"/>	Dual Multi-Axial Headrest hardware -----	OFH3	\$195
<input type="checkbox"/>	Fixed Headrest Hardware Ass'y (Plus 4") -----	OFH4	\$20
<input type="checkbox"/>	Single Multi-Axial Headrest Hardware (Plus 4") -----	OFH5	\$183
<input type="checkbox"/>	Dual Multi-Axial Headrest hardware (Plus 4") -----	OFH6	\$220

HEADREST FOAM PAD OPTIONS

<input type="checkbox"/>	Small Headrest Foam Pad -----	OHP1	No Charge
<input type="checkbox"/>	Small Headrest with Wings Foam Pad -----	OHP2	\$90
<input type="checkbox"/>	Large Headrest Foam Pad -----	OHP3	Standard
<input type="checkbox"/>	Large Headrest with Wings Foam Pad -----	OHP4	\$95
<input type="checkbox"/>	Omit Headrest -----	OHT	less net \$10

MISCELLANEOUS ACCESSORIES



<input type="checkbox"/>	Tilt Locking Collar -----	ATC	\$20
<input type="checkbox"/>	Recline Locking Collar -----	ARC	\$20
<input type="checkbox"/>	Auto Seat Positioning Strap -----	ARB	\$52
<input type="checkbox"/>	Airline Seat Positioning Strap -----	RRB	\$76
<input type="checkbox"/>	Velcro Release Positioning Strap -----	VRB	\$48
<input type="checkbox"/>	Foam Pad on Arm Guard (pair) -----	FPAG	\$140
<input type="checkbox"/>	Clear Lap Tray 15"-18" -----	LTM	\$255
<input type="checkbox"/>	Clear Lap Tray 19"-22" -----	LTN	\$255
<input type="checkbox"/>	Padded Lap Tray 15"-16" -----	PLT16	No Charge
<input type="checkbox"/>	Padded Lap Tray 17"-18" -----	PLT18	No Charge
<input type="checkbox"/>	Padded Lap Tray 19"-20" -----	PLT20	No Charge
<input type="checkbox"/>	Padded Lap Tray 21"-22" -----	PLT22	No Charge
<input type="checkbox"/>	Padded Lap Tray 23"-24" -----	PLT24	No Charge
<input type="checkbox"/>	Half Padded Lap Tray, left side -----	HPLTL	\$160
<input type="checkbox"/>	Half Padded Lap Tray, right side -----	HPLTR	\$160
<input type="checkbox"/>	Elbow Stop, Padded Lap Tray, left side -----	ESPL	\$55
<input type="checkbox"/>	Elbow Stop, Padded Lap Tray, right side -----	ESPR	\$55
12	<input type="checkbox"/>	Shoulder Lateral Pad (each) -----	SHLP \$90
13	<input type="checkbox"/>	Adjustable Side Pad -----	ASP \$157
	<input type="checkbox"/>	IV Pole, telescopic -----	IVP \$168
	<input type="checkbox"/>	Oxygen Tank Holder -----	OTH \$173
	<input type="checkbox"/>	Rear Anti-Tippers -----	ATR Standard
	<input type="checkbox"/>	Quick Release Axle -----	QRA \$110
	<input type="checkbox"/>	Adductor Pad 10W x 4.5H w/Flush Hooks (pair) -----	APSMFL \$175
	<input type="checkbox"/>	Adductor Pad 10W x 4.5H w/Offset Hooks (pair) -----	APSMLOF \$175
	<input type="checkbox"/>	Adductor Pad 12W x 4.5H w/Flush Hooks (pair) -----	APMEDFL \$175
	<input type="checkbox"/>	Adductor Pad 12W x 4.5H w/Offset Hooks (pair) -----	APMEDOF \$175
	<input type="checkbox"/>	Adductor Pad 14W x 4.5H w/Flush Hooks (pair) -----	APLRGFL \$175
	<input type="checkbox"/>	Adductor Pad 14W x 4.5H w/Offset Hooks (pair) -----	APLRGOF \$175
15	<input type="checkbox"/>	Orion II Front and Rear Transit Hooks -----	HOOK \$250
14	<input type="checkbox"/>	Power Tilt -----	ORIONPT \$1,600

FOOTNOTES:

- 12. Only available with Stroller Bar push option.
- 13. Not available with Flip Back Arms
- 14. Weight Capacity of Orion II with Power Tilt is 250 lbs. Power Tilt can only be placed on Orion II chairs with seat width of 15" to 20"
Power Tilt provides approximately 42 degrees of infinite tilt (3° to 45°)
Includes: 2 Button Hand Control, 27V Plugin Battery Charger, 24V Battery/Control Box, Compact Actuator with 135 mm (5.31") Stroke Length
- 15. Transit Hooks can only be used on Orion II wheelchairs up to 350 lbs weight capacity. Tested in accordance with ANSI/RESNA WC/19

WARRANTY

Future Mobility Healthcare Inc. warrants the Orion II to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user. The batteries required for Power Tilt are warranted for six (6) months while the control box, battery charger, actuator and hand control are warranted for one (1) year.



STF MEASURED	REAR WHEEL	MOUNT LEVEL	CASTER WHEEL	FORK SIZE	HOLE #
13.25	20	SUPER-LOW	4"	X-SMALL	ONLY
14.50	20	UPPER	5"	SMALL	1
15.00	20	MIDDLE	5"	SMALL	2
15.50	22	UPPER	5"	SMALL	3
16.00	24	UPPER	6"	SMALL	3
16.00	22	MIDDLE	6"	SMALL	3
16.00	20	LOWER	5"	SMALL	4
16.50	24	UPPER	7"	SMALL	3
16.50	22	MIDDLE	7"	SMALL	3
16.50	24	MIDDLE	6"	SMALL	4
17.00	24	MIDDLE	7"	SMALL	4
17.00	22	LOWER	8"	LARGE	3
18.00	24	LOWER	8"	LARGE	4
19.00	24	LOWER	8"	LARGE	4

NOTE: ALL SEAT TO FLOOR MEASUREMENTS ARE ROUNDED TO THE NEAREST 1/4" 19" STF USES 1" RISER BLOCKS ON SEATPAN