

Phantom Back™

Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market.

The Phantom Back™ is designed to target active wheelchair users, those that require back support and need maximum freedom for upper body movement.

The Phantom Back™ fits rigid wheelchair frames without push handles and replaces standard back support upholstery. It offers maximum seat-to-back infinite angle adjustment within the range, providing real individual fit. Well designed back support panel and cushion promotes lateral and posterior postural support in pelvic region needed for fatigue reduction, thus enhancing upper body function.

The accessibility of hardware interface, the flexibility to adjust and the light weight makes the Phantom Back™ exceptionally meeting the needs for the active user.



SPECIFICATIONS

Height Adjustment:

Infinite adjustment along wheelchair posts height. Hardware and slotted holes in the back panel provide 3" adjustment range in increments of 1", to overcome cross brace obstacle on the wheelchair posts.

Width Adjustment:

Hardware allows width adjustment range of 3/4", for example 16" wide back would fit wheelchair frames from 15 3/4" to 16 1/2".

Angle Adjustment:

In relation to wheelchair back support posts, hardware accommodates up to 15 degrees anterior and up to 15 degrees posterior infinite angle adjustment within this range.

Seat Depth Adjustment: (measurements taken from front side of wheelchair posts)
Optimal seat depth of 1 1/4" behind wheelchair posts available with standard hardware, when Phantom backrest aligned with wheelchair back support posts.

Maximum seat depth is 2 1/2" behind wheelchair posts, when anterior back angle is 15 degrees, minimum seat depth is 0" behind wheelchair posts, when posterior angle is set to 15 degrees.

For shorter seat depth than 0", in addition to standard hardware, a custom spacer is available to set seat depth up to 2" in front of wheelchair posts.

Wheelchair Posts Interface:

Standard hardware will properly fit 1" tubing diameter. For wheelchairs with 7/8" or 3/4" tubing diameter, insert spacers which are included in the hardware.

WARRANTY

Three (3) years limited on Back and 90 days on cover, 1 year on hardware.

PRODUCT WEIGHT

2.7 lbs. (18" wide)

USER WEIGHT

250 lbs.

SIZES

Widths: 12" - 18"



Future Mobility Healthcare

3223 Orlando Drive, Mississauga, ON, L4V 1C5, 905-671-1661, 888-737-4011, fax 866-670-3377 www.future-mobility.com
Future Mobility Healthcare Inc. Orion and Prism are trademarks of Future Mobility Healthcare Inc.