



Non-Powered Adjustable Zoned  
Pressure Reducing Mattress Overlay



The StarMatt Overlay is a four section, therapeutic, air adjustable mattress overlay featuring 180 cells per section for a total of 720 air cells throughout the entire system. As the individual becomes immersed into the StarMatt, they benefit from a high resolution, low pressure, low shear support surface that conforms to the unique contours of the body, providing an ideal environment for the healing of up to and including a Stage IV pressure ulcer. The “**GOLD**” Standard in therapeutic mattress overlays!



A vapor permeable, quilted top cover is included with each StarMatt System. The cover encloses not only the StarMatt, but the hospital or bed mattress as well, protecting the bed mattress from fluids and helping to prevent sliding of the StarMatt during transfers in and out of bed.

A heavy duty carry bag is included to carry the entire StarMatt System, providing an easy method for transport.



StarMatt®  
Features

- Light weight, approximately 32 pounds or 14.5 kilograms per total system weight.
- Non-powered, no noisy pumps or powered components.
- Low maintenance!
- 4” cell height for Maximum Immersion and Pressure Distribution.
- 24 Month Manufacturers Limited Warranty.
- Designed from high quality, flame resistant, latex free neoprene rubber.
- Available in Bariatric sizing up to a width of 60”.
- No weight limit.

**Future Mobility Healthcare**

Future Mobility Healthcare Inc., 3223 Orlando Drive, Mississauga, ON L4V 1C5  
Tel: 1-888-737-4011, Fax: 1-866-670-3377



The StarMatt Replacement System is an alternative to the overlay system. Utilizing a 7" thick foam mattress, a cut-out or well is created to receive the four StarMatt sections. By this design, 3" foam side bolsters are created to ease transfers and aids in patient safety. The design also features soft, flexible air cells throughout the entire foot section of the bed to ensure heel protection. Custom configurations are also available for Bariatric beds as well as King and Queen size mattresses.

STARMATT ®	Description	Suggested Retail
SM 3575	Complete overlay mattress system (34.5" wide)	\$ 3,500
SM 3975	Complete Bariatric overlay mattress system (39" wide)	\$ 3,900
SM 4275	Complete Bariatric overlay mattress system (42" wide)	\$ 4,200
SM 4875	Complete Bariatric overlay mattress system (48" wide)	\$ 4,800
SMR-1	Complete replacement mattress system: includes 1 bed section and 3 foam sections	\$ 3,000
SMR-2	Complete replacement mattress system: includes 2 bed sections and 2 foam sections	\$ 3,200
SMR-3	Complete replacement mattress system: includes 3 bed sections and 1 foam sections	\$ 3,500
SMR-4	Complete replacement mattress system: includes 4 bed sections	\$ 3,750
SM 3520	Single mattress section for overlay (33.5" x 19.25")	\$ 870
SM 3521	Single mattress section for replacement (28" x 19.25")	\$ 870
SM CV01	Quilted mattress cover	\$ 225
SM DFL-1	Carrying bag	\$ 125
SM F1	Single foam section for overlay (33.5" x 19.25")	\$ 250
SM F2	Single foam section for replacement (28" x 19.25")	\$ 250

- Each overlay system includes 4 bed sections, a quilted mattress cover, hand pump, repair kit and carrying bag
- Each replacement system includes a hand pump and repair kit