

# STELLATO

*Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market.*



The Stellato wheelchair is an innovative, low maintenance, performance wheelchair. The Stellato comes with many standard features, including a side storage pocket for personal use that is easily accessible for the user. Developed for comfort and affordability, the Stellato wheelchair is sure to become the model of choice.

**Future M****obility Healthcare**



### FEATURES:

- Lightweight folding frame with swingaway front end
- Elliptical heavy duty crossbrace tubing
- Seat widths 14"-20"
- Seat depths 15"-20"
- Seat to floor height 13 1/4"- 19"
- Rear wheel options 20", 22" or 24"
- Handrim options (Stainless Steel)
- Caster options 4", 5", 7" or 8"
- Camber adjustment caster socket/precision bearings
- Back cane tall- 17-20 Top of seat rail to top edge of back sling
- Back cane short - 14-17 Top of seat rail to top edge of back sling
- Seat upholstery
- Slip on back upholstery
- Angle adjustable push canes
- Aluminum foot plates
- Heel loops
- Wheel lock options—push or pull to lock
- Side storage pocket
- Choice of two colours—Black or Silver



### WARRANTY

Frame is warranted for the lifetime of the original purchase or user.

### PRODUCT WEIGHT

31 lbs.

### USER WEIGHT

250 lbs.

**Future Mobility Healthcare**