

Order Form/Price List

(effective August 18, 2008)

Account Information:

Date: _____

Account #: _____
 Name: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____
 Comment: _____

PO# _____
 Ship To: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____

1 FRAME TYPE

Folding Frame w/Swingaway Front End STELLATO \$1,350

FRAME FINISH

Black FBK No Charge
 Silver FSR No Charge

2 SIDE FRAME MATRIX

SEAT DEPTH						
WIDTH	15	16	17	18	19	20
14	SMALL	SMALL	SMALL	SMALL		
15	SMALL	SMALL	SMALL	SMALL		
16	SMALL	SMALL	SMALL	SM/ LG	LARGE	LARGE
17	SMALL	SMALL	SMALL	SM/ LG	LARGE	LARGE
18	SMALL	SMALL	SMALL	SM/ LG	LARGE	LARGE
19				LARGE	LARGE	LARGE
20				LARGE	LARGE	LARGE

SEAT WIDTH

14" Seat Width SW14 No Charge
 15" Seat Width SW15 \$125
 16" Seat Width SW16 No Charge
 17" Seat Width SW17 No Charge
 18" Seat Width SW18 No Charge
 19" Seat Width SW19 \$125
 20" Seat Width SW20 No Charge

SEAT DEPTH

⁵ 15" Seat Depth SD15 \$125
 16" Seat Depth SD16 No Charge
 17" Seat Depth SD17 \$125
 18" Seat Depth SD18 No Charge
 19" Seat Depth SD19 \$125
 20" Seat Depth SD20 No Charge

3, 6 SEAT-TO-FLOOR HEIGHT - small frame (Hemi)

⁴ 13.25" (20" rear & 4" caster) STFS01 No Charge
 14.5" (20" rear & 5" caster) STFS03 No Charge
 15" (20" rear & 5" caster) STFS04 No Charge
 15.5" (20" rear & 5" caster) STFS05 No Charge
 16" (22" rear & 5" caster) STFS06 No Charge
 16" (22" rear & 6" caster) STFS11 No Charge
 16.5" (22" rear & 7" caster) STFS07 No Charge
 16.5" (22" rear & 6" caster) STFS12 No Charge
 17" (22" rear & 7" caster) STFS08 No Charge
 17.5" (24" rear & 8" caster) STFS09 No Charge
 18" (24" rear & 8" caster) STFS10 No Charge

Basic Chair Includes:

Lightweight Folding Frame w/Swingaway Front End
 Elliptical Heavyduty Crossbrace Tubing
 Seat Width (14"-20")
 Seat Depth (15"-20")
 Seat-To-Floor Height (13.25"-19.0")
 Rear Wheel Options (20", 22" or 24") Airless or Pneumatic
 Handrim Options (Stainless Steel)
 Caster Options (4", 5", 6", 7" or 8")
 Camber Adjustable Caster Socket / Precesion Bearings
 Seat Upholstery
 Slip On Back Upholstery
 Angle Adjustable Push Canes
 Aluminum Flip-Up Footplate
 Heel Loops
 Wheel Locks (Push or Pull to Lock)
 Side Storage Pocket
 Choice of Two Colors (Black or Silver Powdercoat Finish)

FOOTNOTES:

1. Seat angle is set at 3 degrees. See service manual for details
2. For 18" Seat Depth (16"-18"W) both large and small frame available
3. Refer to Side Frame Matrix for SM/LG frame determination
4. Only available in 16" or 18" seat depth
5. Not available with 24" and 8" combination
6. low push canes 16" -19"

13 SEAT-TO-FLOOR HEIGHT - large frame

<input type="checkbox"/> 15.5" (20" rear & 5" caster)	STFL01	No Charge
<input type="checkbox"/> 16" (22" rear & 5" caster)	STFL02	No Charge
<input type="checkbox"/> 16.5" (22" rear & 5" caster)	STFL03	No Charge
<input type="checkbox"/> 16.5" (22" rear & 6" caster)	STFL10	No Charge
<input type="checkbox"/> 17" (22" rear & 7" caster)	STFL04	No Charge
<input type="checkbox"/> 17" (24" rear & 6" caster)	STFL11	No Charge
<input type="checkbox"/> 17.5" (24" rear & 7" caster)	STFL05	No Charge
<input type="checkbox"/> 17.5" (24" rear & 6" caster)	STFL12	No Charge
<input type="checkbox"/> 18" (24" rear & 7"caster)	STFL06	No Charge
<input type="checkbox"/> 18" (24" rear & 6"caster)	STFL13	No Charge
<input type="checkbox"/> 18.5" (24" rear & 7" caster)	STFL07	No Charge
<input type="checkbox"/> 18.5" (24" rear & 8" caster)	STFL08	No Charge
<input type="checkbox"/> 19" (24" rear & 8" caster)	STFL09	No Charge

FOOTNOTES:

- 7. Foam pads are recommended if using lap tray.
- 8. Back angle adjusts -7 degrees to 21 degrees in 7 degree increments
- 9. Not available with composite footplate with 6", 7" or 8" casters.
- 10. Not available with angle adjustable footplate with 6", 7" or 8" casters.
- 11. Not available with 6", 7" or 8" caster when footboard is facing back on 90 degree swingaway.
- 12. Please refer to Backs order form for stock,non stock and custom sizes
- 13. high push canes 18"-21"

WHEEL OPTIONS

<input type="checkbox"/> Airless	WHLA	Standard
<input type="checkbox"/> Pneumatic	WHLP	No Charge

HANDRIM OPTIONS

<input type="checkbox"/> Stainless Steel	SS	Standard
<input type="checkbox"/> Plastic Coated	SPC	\$110
<input type="checkbox"/> Omit Handrims	OHR	less net \$15

ARM TYPE

7 <input type="checkbox"/> Full Length Adj. Height Flipback (foam pad)	ARMCF05	\$242
<input type="checkbox"/> Full Length Adj. Height Flipback (gel pad)	ARMCG05	\$298
<input type="checkbox"/> Full Length Fixed Height Flipback (vinyl pad)	ARMFH	\$195

BACK TYPE

8 <input type="checkbox"/> Angle Adj 10 deg Bend w/Push Handle (low)	SBCN110	Standard
8 <input type="checkbox"/> Angle Adj 10 deg Bend w/Push Handle (high)	SBCN111	\$90

BACK UPHOLSTERY

<input type="checkbox"/> Slip-on Padded Back Upholstery	BUSB	Standard
<input type="checkbox"/> Adjustable Tension Back	TB	\$175
<input type="checkbox"/> Prism Basic Back on chair Stock -----	SBBOCS	\$500
<input type="checkbox"/> Prism Basic Back on chair Non Stock -----	SBBOCNS	\$525
<input type="checkbox"/> Prism Basic Back on chair Cus -----	SBBCUS	\$825
<input type="checkbox"/> Prism Ultra Back on chair Stock -----	SUBOCS	\$600
<input type="checkbox"/> Prism Ultra Back on chair Non Stock -----	SUBOCNS	\$625
<input type="checkbox"/> Prism Ultra Back on chair Cus -----	SUBCUS	\$1,010
<input type="checkbox"/> Prism Air Back on chair Stock -----	SABOCS	\$650
<input type="checkbox"/> Prism Air Back on chair Non Stock -----	SABOCNS	\$675
<input type="checkbox"/> Prism Air Back on chair Cus -----	SABCUS	\$1,075
<input type="checkbox"/> Prism TrueFitt Back on chair Stock -----	STFBOCS	\$750
<input type="checkbox"/> Prism TrueFitt Back on chair Non Stock -----	STFBOCNS	\$775
<input type="checkbox"/> Prism TrueFitt Back on chair Cus -----	STFBCUS	\$1,025
<input type="checkbox"/> Omit Back Upholstery	OBU	less net \$20
12 <input type="checkbox"/> Specify Back Size W____x H____		

SEAT OPTIONS

<input type="checkbox"/> Black Nylon Seat Upholstery	SUN	Standard
<input type="checkbox"/> Omit Seat Upholstery	OSU	less net \$20
<input type="checkbox"/> Prism Ideal Cushion - Stock Sizes	ICSCS	\$175
<input type="checkbox"/> Prism Ideal Cushion - Non-Stock Sizes	ICSCNS	\$200
<input type="checkbox"/> Prism Ideal Cushion - Custom Sizes	ICCUSOC	\$275
<input type="checkbox"/> Prism Elite Cushion - Stock Sizes	ECOC	\$250
<input type="checkbox"/> Prism Elite Cushion - Non-Stock Sizes	ECOCNS	\$275
<input type="checkbox"/> Prism Elite Cushion - Custom Sizes	ECCUSOC	\$425
<input type="checkbox"/> Prism Supreme Cushion - Stock Sizes	SCOC	\$375
<input type="checkbox"/> Prism Supreme Cushion - Non-Stock Sizes	SCOCNS	\$400
<input type="checkbox"/> Prism Supreme Cushion - Custom Sizes	SCCUSOC	\$685
<input type="checkbox"/> Prism Supreme II Plus Cushion Stock on cha	SCPOC	\$375
<input type="checkbox"/> Prism Supreme II Plus Cushion Non-Stock oi	SCPOC	\$400
<input type="checkbox"/> Prism Supreme II Plus Cushion Custom on c	SCPCUS	\$685
<input type="checkbox"/> Prism Supreme Air Cushion - Stock Sizes	SACOC	\$500
<input type="checkbox"/> Prism Supreme Air Cushion - Non-Stock Sizes	SACOCNS	\$575
<input type="checkbox"/> Prism Supreme Drop Base 1" Stock Sizes	SC1DB	\$755
<input type="checkbox"/> Prism Supreme Drop Base 2" Stock Sizes	SC2DB	\$755

- Prism Supreme Drop Base 1" Non Stock SC1DBNS \$810
- Prism Supreme Drop Base 2" Non Stock SC2DBNS \$810
- Prism Drop Base 1" Stock Sizes DB1 \$180
- Prism Drop Base 2" Stock Sizes DB2 \$180
- Prism Drop Base 1" Non Stock Sizes DB1NS \$210
- Prism Drop Base 2" Non Stock Sizes DB2NS \$210

FRONT RIGGING

- 70 Degree Swingaway Pinstyle SFRP70 \$132
- 80 Degree Swingaway Pinstyle SFRP80 \$135
- 60 Degree Swingaway Upgraded SF60 No Charge
- 70 Degree Swingaway Upgraded SF70 No Charge
- ^{9,10} 90 Degree Swingaway Upgraded SF90 \$230
- Elevating Legrest (pair) w/adj. Height Hanger ELR \$550

FOOTPLATE TYPE

- ⁹ Composite Flip-Up Footplates FPC No Charge
- Aluminum Flip-Up Footplates FPA Standard
- X-Large Aluminum Flip-Up Footplates XFPA \$90
- ¹⁰ Adjustable Angle Flip-Up Footplates FFAA \$120
- ¹¹ One Piece Flip-Up Footboard FPOP \$160

FOOTREST/LEGREST ACCESSORIES

- Padded Calf Strap PCS \$70
- Foam Guards on Elevating Legrest (pair) FG \$70
- Heel Loops HL Standard
- Omit Heel Loops OHL less net \$7

WHEEL LOCKS

- Push-to-Lock Wheel Lock WLSPL Standard
- Pull-to-Lock Wheel Lock WLSPU No Charge
- Wheel Lock Extension WLEXT No Charge

MISCELLANEOUS ACCESSORIES

- Anti-tippers ATS No Charge
- Auto-Seat Positioning Straps ARB \$52
- Velcro Release Positioning Strap VRB \$48
- Quick Release Axles QRA \$110
- Clear Lap Tray 15"-18" LTM \$255
- Clear Lap Tray 19"-22" LTL \$255
- Padded Lap Tray 15"-16" SPLT16 \$230
- Padded Lap Tray 17"-18" SPLT18 \$230
- Padded Lap Tray 19"-20" SPLT20 \$230
- Padded Lap Tray 21"-22" SPLT22 \$230
- Padded Lap Tray 23"-24" SPLT24 \$230
- Half Padded Lap Tray, left side SHPLTL \$160
- Half Padded Lap Tray, right side SHPLTR \$160
- IV Pole, telescopic SIVP \$168
- Oxygen Tank Holder SOTH \$175

SEAT-TO-FLOOR HEIGHT MATRIX

STF	REAR	REAR	CASTER	FORK	FORK	MOUNT
	WHEEL	HOLE	WHEEL	SIZE	HOLE	POS.
SMALL FRAME						
13.25	20	1	4	X-SM	1	1
14.50	20	2	5	SHORT	1	1
15.00	20	3	5	SHORT	2	1
15.50	20	3	5	SHORT	3	1
16.00	22	2	5	LONG	2	1
16.00	22	3	6	SHORT	3	1
16.50	22	3	7	SHORT	3	1
16.50	22	4	6	LONG	2	1
17.00	22	4	7	LONG	2	1
17.50	24	3	8	LONG	2	1
18.00	24	4	8	LONG	3	1
LARGE FRAME						
15.50	20	3	5	SHORT	1	1
16.00	22	2	5	SHORT	2	1
16.50	22	3	5	SHORT	3	1
16.50	22	4	6	LONG	2	1
17.00	22	4	7	LONG	2	1
17.00	24	3	6	SHORT	3	1
17.50	24	3	7	SHORT	3	1
17.50	24	4	6	LONG	2	1
18.00	24	4	7	LONG	2	1
18.00	24	4	6	LONG	3	1
18.50	24	4	7	LONG	3	1
18.50	24	4	8	LONG	2	1
19.00	24	5	8	LONG	3	1