

Order Form/Price List

(Effective August 6, 2010)

Account Information:

Date: _____

Account #: _____
 Name: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____
 Comment: _____

PO# _____
 Ship To: _____
 Address: _____
 City: _____
 Province: _____
 Postal: _____
 Phone: _____
 Fax: _____

1 FRAME TYPE

- Folding Frame w/Swingaway Front End STELLATO \$1,350

FRAME FINISH

- Black FBK No Charge
- Silver FSR No Charge
- Custom Colours FCUS NET \$150

2 SIDE FRAME MATRIX

SEAT DEPTH						
WIDTH	15	16	17	18	19	20
14	SMALL	SMALL	SMALL	SMALL		
15	SMALL	SMALL	SMALL	SMALL		
16	SMALL	SMALL	SMALL	SM/ LG	LARGE	LARGE
17	SMALL	SMALL	SMALL	SM/ LG	LARGE	LARGE
18	SMALL	SMALL	SMALL	SM/ LG	LARGE	LARGE
19				LARGE	LARGE	LARGE
20				LARGE	LARGE	LARGE

SEAT WIDTH

- 14" Seat Width SW14 No Charge
- 15" Seat Width SW15 \$125
- 16" Seat Width SW16 No Charge
- 17" Seat Width SW17 No Charge
- 18" Seat Width SW18 No Charge
- 19" Seat Width SW19 \$125
- 20" Seat Width SW20 No Charge

SEAT DEPTH

- 15" Seat Depth SD15 \$125
- 16" Seat Depth SD16 No Charge
- 17" Seat Depth SD17 \$125
- 18" Seat Depth SD18 No Charge
- 19" Seat Depth SD19 \$125
- 20" Seat Depth SD20 No Charge

3, 6 SEAT-TO-FLOOR HEIGHT - small frame (Hemi)

- 13.25" (20" rear & 4" caster) STFS01 No Charge
- 14.5" (20" rear & 5" caster) STFS03 No Charge
- 15" (20" rear & 5" caster) STFS04 No Charge
- 15.5" (20" rear & 5" caster) STFS05 No Charge
- 16" (22" rear & 5" caster) STFS06 No Charge
- 16" (22" rear & 6" caster) STFS11 No Charge
- 16.5" (22" rear & 7" caster) STFS07 No Charge
- 16.5" (22" rear & 6" caster) STFS12 No Charge
- 17" (22" rear & 7" caster) STFS08 No Charge
- 17.5" (24" rear & 8" caster) STFS09 No Charge
- 18" (24" rear & 8" caster) STFS10 No Charge

Basic Chair Includes:

Lightweight Folding Frame w/Swingaway Front End
 Elliptical Heavyduty Crossbrace Tubing
 Seat Width (14"-20")
 Seat Depth (15"-20")
 Seat-To-Floor Height (13.25"-19.0")
 Rear Wheel Options (20", 22" or 24") Airless or Pneumatic
 Handrim Options (Stainless Steel)
 Caster Options (4", 5", 6", 7" or 8")
 Camber Adjustable Caster Socket / Precision Bearings
 Seat Upholstery
 Slip On Back Upholstery
 Angle Adjustable Push Canes
 Aluminum Flip-Up Footplate
 Heel Loops
 Wheel Locks (Push or Pull to Lock)
 Side Storage Pocket
 Choice of Two Colors (Black or Silver Powdercoat Finish)

FOOTNOTES:

1. Seat angle is set at 3 degrees. See service manual for details
2. For 18" Seat Depth (16"-18"W) both large and small frame available
3. Refer to Side Frame Matrix for SML/LG frame determination
4. Only available in 16" or 18" seat depth
5. Not available with 24" and 8" combination
6. low push canes 16" -19"

13 SEAT-TO-FLOOR HEIGHT - large frame

<input type="checkbox"/>	15.5" (20" rear & 5" caster)	STFL01	No Charge
<input type="checkbox"/>	16" (22" rear & 5" caster)	STFL02	No Charge
<input type="checkbox"/>	16.5" (22" rear & 5" caster)	STFL03	No Charge
<input type="checkbox"/>	16.5" (22" rear & 6" caster)	STFL10	No Charge
<input type="checkbox"/>	17" (22" rear & 7" caster)	STFL04	No Charge
<input type="checkbox"/>	17" (24" rear & 6" caster)	STFL11	No Charge
<input type="checkbox"/>	17.5" (24" rear & 7" caster)	STFL05	No Charge
<input type="checkbox"/>	17.5" (24" rear & 6" caster)	STFL12	No Charge
<input type="checkbox"/>	18" (24" rear & 7"caster)	STFL06	No Charge
<input type="checkbox"/>	18" (24" rear & 6"caster)	STFL13	No Charge
<input type="checkbox"/>	18.5" (24" rear & 7" caster)	STFL07	No Charge
<input type="checkbox"/>	18.5" (24" rear & 8" caster)	STFL08	No Charge
<input type="checkbox"/>	19" (24" rear & 8" caster)	STFL09	No Charge

WHEEL OPTIONS

<input type="checkbox"/>	Airless	WHLA	Standard
<input type="checkbox"/>	Pneumatic	WHLP	No Charge
<input type="checkbox"/>	Pneumatic Spoke Wheels	WHLS	\$175

HANDRIM OPTIONS

<input type="checkbox"/>	Stainless Steel	SS	Standard
<input type="checkbox"/>	Plastic Coated	SPC	\$110
<input type="checkbox"/>	Omit Handrims	OHR	less net \$15

ARM TYPE

7	<input type="checkbox"/>	Full Length Adj. Height Flipback (foam pad)	ARMCFO5	\$242
	<input type="checkbox"/>	Full Length Adj. Height Flipback (gel pad)	ARMCGO5	\$298
	<input type="checkbox"/>	Full Length Fixed Height Flipback (vinyl pad)	ARMFH	\$195

BACK TYPE

8	<input type="checkbox"/>	Angle Adj 10 deg Bend w/Push Handle (low)	SBCN110	Standard
8	<input type="checkbox"/>	Angle Adj 10 deg Bend w/Push Handle (high)	SBCN111	\$90

BACK UPHOLSTERY

	<input type="checkbox"/>	Slip-on Padded Back Upholstery	BUSB	Standard
	<input type="checkbox"/>	Adjustable Tension Back	TB	\$175
	<input type="checkbox"/>	Prism Basic Back on chair Stock	SBBOC	\$500
	<input type="checkbox"/>	Prism Basic Back on chair Non Stock	SBBOC	\$525
	<input type="checkbox"/>	Prism Basic Back on chair Cus	SBBCUS	\$825
	<input type="checkbox"/>	Prism Ultra Back on chair Stock	SUBOC	\$600
	<input type="checkbox"/>	Prism Ultra Back on chair Non Stock	SUBOC	\$625
	<input type="checkbox"/>	Prism Ultra Back on chair Cus	SUBCUS	\$1,010
	<input type="checkbox"/>	Prism Air Back on chair Stock	SABOC	\$650
	<input type="checkbox"/>	Prism Air Back on chair Non Stock	SABOC	\$675
	<input type="checkbox"/>	Prism Air Back on chair Cus	SABCUS	\$1,075
	<input type="checkbox"/>	Prism TrueFitt Back on chair Stock	STFBOC	\$750
	<input type="checkbox"/>	Prism TrueFitt Back on chair Non Stock	STFBOC	\$775
	<input type="checkbox"/>	Prism TrueFitt Back on chair Cus	STFBCUS	\$1,025
	<input type="checkbox"/>	Prism Max Ultra Back on chair Stock	SMUBOC	\$650
	<input type="checkbox"/>	Prism Max Ultra Back on chair Non Stock	SMUBOC	\$675
	<input type="checkbox"/>	Prism Max Ultra Back on chair Cus	SMUBCUS	\$1,060
	<input type="checkbox"/>	Prism Max Air Back on chair Stock	SMABOC	\$700
	<input type="checkbox"/>	Prism Max Air Back on chair Non Stock	SMABOC	\$725
	<input type="checkbox"/>	Prism Max Air Back on chair Cus	SMABCUS	\$1,125
	<input type="checkbox"/>	Prism Max TrueFitt Back on chair Stock	SMTFBOC	\$800
	<input type="checkbox"/>	Prism Max TrueFitt Back on chair Non Stock	SMTFBOC	\$825
	<input type="checkbox"/>	Prism Max TrueFitt Back on chair Cus	SMTFBCUS	\$1,075
14	<input type="checkbox"/>	Lunar Back Low No Lateral on chair	SLBLNOC	\$500
	<input type="checkbox"/>	Lunar Back Low Basic Lateral on chair	SLBLBOC	\$575
15	<input type="checkbox"/>	Lunar Back Mid No Lateral on chair	SLBMNOC	\$590
	<input type="checkbox"/>	Lunar Back Mid Basic Lateral on chair	SLBMBOC	\$650
	<input type="checkbox"/>	Lunar Back Mid Ultra Lateral on chair	SLBMUOC	\$710
16	<input type="checkbox"/>	Lunar Back Full No Lateral on chair	SLBFNOC	\$610
	<input type="checkbox"/>	Lunar Back Full Basic Lateral on chair	SLBFBOC	\$670
	<input type="checkbox"/>	Lunar Back Full Ultra Lateral on chair	SLBFUOC	\$730
	<input type="checkbox"/>	Lunar Back Custom on chair	SLBCUSON	\$1,050
	<input type="checkbox"/>	Omit Back Upholstery	OBU	less net \$20
12	<input type="checkbox"/>	Specify Back Size W____x H____		

FOOTNOTES:

- 7. Foam pads are recommended if using lap tray.
- 8. Back angle adjusts -7 degrees to 21 degrees in 7 degree increments
- 9. Not available with composite footplate with 6", 7" or 8" casters.
- 10. Not available with angle adjustable footplate with 6", 7" or 8" casters.
- 11. Not available with 6", 7" or 8" caster footboard is facing back on 90 degree swingaway.
- 12. Please refer to Backs order form for stock, non stock and custom sizes
- 13. high push canes 18"-21"
- * 14. Sizes 14W, 16W, 18W, 20W x 10H or 12H
- * 15. Sizes 14W, 16W, 18W, 20W x 14H, 16H or 18H
- * 16. Sizes 16W, 18W, 20W, 22W x 20H, 22H or 24H
- * (refer to lunar backs order form for specific sizes)

SEAT OPTIONS

<input type="checkbox"/>	Black Nylon Seat Upholstery	SUN	Standard
<input type="checkbox"/>	Omit Seat Upholstery	OSU	less net \$20
<input type="checkbox"/>	Prism Ideal Cushion - Stock Sizes	SICOC	\$175
<input type="checkbox"/>	Prism Ideal Cushion - Non-Stock Sizes	SICOC	\$200
<input type="checkbox"/>	Prism Ideal Cushion - Custom Sizes	ICCUS	\$275
<input type="checkbox"/>	Prism Elite Cushion - Stock Sizes	ECOC	\$250
<input type="checkbox"/>	Prism Elite Cushion - Non-Stock Sizes	ECOC	\$275
<input type="checkbox"/>	Prism Elite Cushion - Custom Sizes	ECCUS	\$425
<input type="checkbox"/>	Prism Supreme Cushion - Stock Sizes	SCOC	\$375
<input type="checkbox"/>	Prism Supreme Cushion - Non-Stock Sizes	SCOC	\$400
<input type="checkbox"/>	Prism Supreme Cushion - Custom Sizes	SCCUS	\$685
<input type="checkbox"/>	Prism Supreme II Plus Cushion Stock Sizes	SCPOC	\$375
<input type="checkbox"/>	Prism Supreme II Plus Cushion Non-Stock Sizes	SCPOC	\$400
<input type="checkbox"/>	Prism Supreme II Plus Cushion Custom Sizes	SCPCUS	\$685
<input type="checkbox"/>	Prism Supreme Air Cushion - Stock Sizes	SACOC	\$500
<input type="checkbox"/>	Prism Supreme Air Cushion - Non-Stock Sizes	SACOC	\$575
<input type="checkbox"/>	Prism Supreme Air Cushion - Custom Sizes	SACCUS	\$785
<input type="checkbox"/>	Prism Supreme Drop Base Flushmount Stock Sizes	SDB0OC	\$755
<input type="checkbox"/>	Prism Supreme Drop Base 1" Stock Sizes	SDB1OC	\$755
<input type="checkbox"/>	Prism Supreme Drop Base 2" Stock Sizes	SDB2OC	\$755
<input type="checkbox"/>	Prism Supreme Drop Base Flushmount Non Stock	SDB0OC	\$810
<input type="checkbox"/>	Prism Supreme Drop Base 1" Non Stock	SDB1OC	\$810
<input type="checkbox"/>	Prism Supreme Drop Base 2" Non Stock	SDB2OC	\$810
<input type="checkbox"/>	Prism Supreme Drop Base Flushmount Custom Sizes	SDB0CUS	\$935
<input type="checkbox"/>	Prism Supreme Drop Base 1" Custom Sizes	SDB1CUS	\$935
<input type="checkbox"/>	Prism Supreme Drop Base 2" Custom Sizes	SDB2CUS	\$935
<input type="checkbox"/>	Prism Drop Base Flushmount Stock Sizes	DB0OC	\$180
<input type="checkbox"/>	Prism Drop Base 1" Stock Sizes	DB1OC	\$180
<input type="checkbox"/>	Prism Drop Base 2" Stock Sizes	DB2OC	\$180
<input type="checkbox"/>	Prism Drop Base Flushmount Non Stock Sizes	DB0OC	\$210
<input type="checkbox"/>	Prism Drop Base 1" Non Stock Sizes	DB1OC	\$210
<input type="checkbox"/>	Prism Drop Base 2" Non Stock Sizes	DB2OC	\$210
<input type="checkbox"/>	Prism Drop Base Flushmount Custom Sizes	DB0CUS	\$250
<input type="checkbox"/>	Prism Drop Base 1" Custom Sizes	DB1CUS	\$250
<input type="checkbox"/>	Prism Drop Base 2" Custom Sizes	DB2CUS	\$250

FRONT RIGGING

<input type="checkbox"/>	70 Degree Swingaway Pinstyle	SFRP70	\$132
<input type="checkbox"/>	80 Degree Swingaway Pinstyle	SFRP80	\$135
<input type="checkbox"/>	60 Degree Swingaway Upgraded	SF60	No Charge
<input type="checkbox"/>	70 Degree Swingaway Upgraded	SF70	No Charge
^{9,10} <input type="checkbox"/>	90 Degree Swingaway Upgraded	SF90	\$230
<input type="checkbox"/>	Elevating Legrest (pair) w/adj. Height Hanger	ELR	\$550

FOOTPLATE TYPE

⁹ <input type="checkbox"/>	Composite Flip-Up Footplates	FPC	No Charge
<input type="checkbox"/>	Aluminum Flip-Up Footplates	FPA	Standard
<input type="checkbox"/>	X-Large Aluminum Flip-Up Footplates	XFPA	\$90
¹⁰ <input type="checkbox"/>	Composite Adjustable Angle Flip-Up Footplates	FFAAC	\$120
¹⁰ <input type="checkbox"/>	Aluminum Adjustable Angle Flip-Up Footplates	FFAA	\$120
¹⁰ <input type="checkbox"/>	Composite Adjustable Angle Flip-Up Footplates Spring Return	FFAASR	\$150
¹¹ <input type="checkbox"/>	One Piece Flip-Up Footboard	FPOP	\$160

FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/>	Padded Calf Strap	PCS	\$70
<input type="checkbox"/>	Foam Guards on Elevating Legrest (pair)	FG	\$70
<input type="checkbox"/>	Foam Pad on Footplate (pair)	FPF	\$70
<input type="checkbox"/>	Foam Pad on One Piece Footboard	OPFD	\$50
<input type="checkbox"/>	Heel Loops	HL	Standard
<input type="checkbox"/>	Omit Heel Loops	OHL	less net \$7

WHEEL LOCKS

<input type="checkbox"/>	Push-to-Lock Wheel Lock	WLSPL	Standard
<input type="checkbox"/>	Pull-to-Lock Wheel Lock	WLSPU	No Charge
<input type="checkbox"/>	Wheel Lock Extension	WLEXT	No Charge

MISCELLANEOUS ACCESSORIES

<input type="checkbox"/>	Anti-tippers	ATS	No Charge
--------------------------	--------------	-----	-----------

