




*Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market.*

The Prism™ Super Back for Kids is a junior size seat back which provides children excellent pressure relief.

Extras, such as, the positioning option, the ease of attachment and the unique method of adjusting for a wide range of angles, makes this Back the most versatile and truly modular model on the market.

Combine this with Prism's sleek and attractive tailored look for a Back that will appeal to a wide variety of users!





**SPOTLIGHT ON**

- Lightweight, sturdy aluminum shell
- Full 5" range of forward and rear adjustability (Angle adjustment forward 15°, back 20°)
- Attractive cover design
- Cane clamps included
- Replaceable foam
- Lumbar pad and contour pad come as standard
- Quick push button release for easy installation

**WARRANTY**

Three (3) years limited on Back and 90 days on cover, 1 year on hardware.

**PRODUCT WEIGHT**

3 lbs

**USER WEIGHT**

150 lbs.

**Stock Sizes—Same day shipments**

16" w x 16" h	KSB1616
14" w x 14" h	KSB1414
12" w x 12" h	KSB1212
10" w x 12" h	KSB1012

**Non-Stock Sizes—Shipments within 3 business days**

12" w x 16" h	KSB1216	11" w x 12" h	KSB1112
13" w x 16" h	KSB1316	13" w x 12" h	KSB1312
14" w x 16" h	KSB1416	14" w x 12" h	KSB1412
11" w x 14" h	KSB1114	15" w x 12" h	KSB1512
12" w x 14" h	KSB1214	16" w x 12" h	KSB1612
15" w x 14" h	KSB1514		
16" w x 14" h	KSB1614		

*Custom orders will rarely be turned away.....we can accommodate as low as 11" in back height and ship it within 7 days!*

