



## Supreme Air Cushion

The Prism Supreme Air Cushion was designed for users that require moderate postural support and, similar to our original Supreme Cushion, offers excellent pressure distribution.

Using inflatable air inserts allows the user to adjust the seat cushion in four different areas for proper weight distribution and comfort. Using our anatomically designed Supreme Cushion allows us to keep our superior pelvic stability and leg positioning, with air inserts for pressure relief.

Supreme comfort, supreme positioning and supreme pressure relief—that's the Prism Supreme Air Cushion!

### SPOTLIGHT ON

- Waterfall front edge for foot propelling
- Four air inserts for optimum comfort and positioning.
- Core layer of visco-elastic foam to provide comfort and pressure reduction

### Prism Supreme Air Cushion

Stock Sizes—same day shipments

14" x 14"	SCA1414	14" x 16"	SCA1416
16" x 16"	SCA1616	16" x 18"	SCA1618
18" x 16"	SCA1816	18" x 18"	SCA1818
18" x 20"	SCA1820		

Non-stocked Sizes—shipments within 3 business days

15" x 16"	SCA1516	15" x 18"	SCA1518
16" x 20"	SCA1620	17" x 16"	SCA1716
17" x 18"	SCA1718	17" x 20"	SCA1720
19" x 16"	SCA1916	19" x 18"	SCA1918
19" x 20"	SCA1920	20" x 16"	SCA2016
20" x 18"	SCA2018	20" x 20"	SCA2020

### WARRANTY


Two (2) years limited on Cushion and 90 days on cover.

### PRODUCT WEIGHT

2.5 lbs.

### USER WEIGHT

250 lbs.



### CUSHION THICKNESS

A—Height at adductor = 5"

B—Height at abductor = 5"

C—Height at rear IT = 3 1/4"

D—Leg trough = 3 1/4"

E—Height at rear trochanter = 4 1/4"

Air inserts = adjustable areas



