

**New**

**PRISM Seating & Positioning**

**Supreme II Plus Cushion**




The Prism Supreme II Plus Cushion was designed for users requiring high positioning support and provides excellent pressure distribution relief.

The anatomically designed firm polyurethane layer with its aggressive leg troughs and abductor contouring provides superior pelvic stability and leg positioning. A base layer of pressure relieving visco-elastic foam under the IT's allows for excellent immersion for bony prominences.

Supreme comfort, supreme positioning and supreme pressure relief—that's the Prism Supreme II Plus!

**SPOTLIGHT ON**

- High density, contour cut layered foam for a consistent, long lasting product
- Waterfall front edge for foot propelling comfort
- Cane cut in rear to allow for correct placement in wheelchairs and proper seat depth
- Radius cut to match shape of back upholstery
- Layer of visco foam along the bottom for further pressure relief under IT's
- Supreme pressure relieving
- Standard water resistant outer cover



**CUSHION THICKNESS**

- A—Height at adductor = 5"
- B—Height at abductor = 5"
- C—Height at rear IT = 3 1/4"
- D—Leg trough = 3 1/4"
- E— Height at rear trochanter = 4 1/4"

**WARRANTY**

Two (2) years limited on Cushion and 90 days on cover.

**PRODUCT WEIGHT**

3.65lbs

**USER WEIGHT**

350lbs



**Prism Supreme Cushion**

*Stock Sizes—same day shipments*

14" x 14"	SCP1414	14" x 16"	SCP1416
16" x 16"	SCP1616	16" x 18"	SCP1618
18" x 16"	SCP1816	18" x 18"	SCP1818
18" x 20"	SCP1820		

*Non-stocked Sizes—shipments within 3 business days*

15" x 16"	SCP1516	15" x 18"	SCP1518
16" x 20"	SCP1620	17" x 16"	SCP1716
17" x 18"	SCP1718	17" x 20"	SCP1720
19" x 16"	SCP1916	19" x 18"	SCP1918
19" x 20"	SCP1920	20" x 16"	SCP2016
20" x 18"	SCP2018	20" x 20"	SCP2020