

## **TrueFitt Back™**

*Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market.*

The Prism™ TrueFitt Back was designed for users with a variety of positioning needs such as Kyphosis, Lordosis, Scoliosis or protruding gibbous deformities.

The ability to customize the Back to each individual's needs, plus the exclusive Prism hardware, allowing for ease of adjustment and installation, contributes to the unique and personal "true fit" in our TrueFitt Back!





### SPOTLIGHT ON

- 1/2" of comfort foam layered on a 1/2" of support foam covering "ice cube" pattern of standard 2" x 2" foam blocks
- Extra foam covering prevents any exposure of shell to user
- Remove, cut or carve individual blocks or entire panels to customize back to fit users' deformities
- "mix and match" foam cubes by using 2" polyurethane foam and ordering our optional 3" polyurethane or 2" visco foam cubes!
- Seat angle adjustability from 82 degrees to 114 degrees optimizing users' positioning allowing for visual orientation and balance
- Lightweight, sturdy aluminum shell

### WARRANTY


Three (3) years limited on Back and 90 days on cover, 1 year on hardware.

### PRODUCT WEIGHT

5.2 lbs. (18" x 16")

### USER WEIGHT

250 lbs. (350 lbs. for bariatric sizes)



### Stock Sizes—Same day shipments

16" x 16"	TFB1616
18" x 16"	TFB1816
20" x 16"	TFB2016

### Non-Stock Sizes—Shipments within 3 business days

17" x 16"	TFB1716
16" x 18"	TFB1618
18" x 18"	TFB1818
16" x 20"	TFB1620
17" x 20"	TFB1720
18" x 20"	TFB1820
20" x 18"	TFB2018
20" x 20"	TFB2020

*Custom orders will rarely be turned away.....we can accommodate as low as 11" in back height and ship it within 7 days!*

