

HEAVY DUTY



Future Mobility Healthcare Inc., manufactures and markets products that increase the freedom and safety of mobility challenged individuals. Reliable for quality and service, Future Mobility Healthcare Inc. provides high quality and innovative mobility products to the home healthcare market.



Orion // 500 is the most innovative and affordable tilt-in-space wheelchair on the market today for people weighting from 250 lbs to 500 lbs. Developed for those individuals requiring adjustability to redistribute their body weight for optimal positioning purposes, the Orion // 500 assists in allowing the user to function independently and to become more involved in their lifestyles.

With a true 25 degrees of tilt, standard on-chair adjustable seating system, incredible range of seat-to-floor heights and creative streamlined design.

The Orion // 500 ... dynamic, progressive, sure to become the most popular heavy duty tilt chair under the stars!



SPOTLIGHT ON:

- Stroller handles
- Standard seating system consisting of: two-way contoured, padded headrest, Orion Back with 1" depth adjustability and Prism Ideal Cushion
- Built in adjustable stop for tilt
- "anti-rattle" anti-tippers
- Quick adjustable contoured armrest pads
- Quick removable arm
- Multi-adjustable axle plate
- Spring return push-to-lock wheel lock
- Unique push button release on optional front rigging
- Innovative elliptical tubing
- Clean, sleek appearance

FEATURES:

- 25 degrees of infinite tilt
- Fixed Back Angle Adjustment: +7 to -21 degrees
- Seat widths: 24"-32" (500 lbs)
- Seat depths: 15"-22"
- Seat-to-floor heights: 15-1/2" - 18"
- Arm heights: 9-1/2" - 14" (adjustable)
- Front rigging: 70 degree standard
- Fixed Back
- 12", 20", 22" or 24" rear wheels
- 5", 6" or 8" casters
- Tall push cane heights
- Aluminum seat pan

SPECIFICATIONS

Tilt Range: 3-28 degrees
 Fixed Back Angle Adjustment: +7 to - 21 degrees

OVERALL DIMENSIONS

(26" x 18" tilt, 16 1/2" STF)
 Chair width: 35"
 Chair length: 49" (with footboard)
 Chair height: 39"
 Chair weight: 70 lbs (not including cushion, back, arm rests or front riggings)

WARRANTY

Frame is warranted for the lifetime of the original purchaser or user.

