


STARLOCK™

Patent Pending

First in the development of complete “CELL LOCK” technology, Star Cushion Products offers you the StarLock Cushion, our premier wheelchair cushion offering you the ultimate in aggressive positioning, unparalleled stability and the highest level of therapeutic prevention and treatment of pressure ulcers you’ve come to expect from a Star Cushion product.



StarLock Features:

- Available in Standard (4”), Mid (3”), and Short (2”) Contour cell heights to better accommodate your seating and positioning needs.
- Variable air flow control provides versatility when addressing hard to fit clients.
- “LOCK TECHNOLOGY” allows you to easily address asymmetrical seating issues.
- “LOCK TECHNOLOGY” allows you to isolate individual air cells, and customize with “no contact” areas without the added expense of an actual custom cushion.
- “LOCK TECHNOLOGY” provides added protection in those cases of accidental cell damage. When locked, only air from damaged cell(s) lose air, not the entire cushion or compartment.
- The ability to isolate air flow between cells enhances stability and assists in transfers and overall posture or trunk control.
- Maintains skin integrity by providing a therapeutic, low pressure, low shear, low friction seating environment.
- Cost effective seating option for a wide range of seating and positioning needs.
- Light weight, weighing approximately 3 pounds/1.36 kilograms.
- Available in Pediatric and Adult Sizes
- 24 Month Manufacturers Limited Warranty
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene rubber.



Leaving a Lasting Impression

The StarLock Cushion combines our traditional/proven air adjustable, multi-cellular technology with our unique “patent pending” air flow control feature that provides aggressive positioning and stability. This innovative system consists of multiple air locking mechanisms called “Lockpockets” located at the base of the cushion. As air volume in the “Lockpockets” is increased through pumping of the attached bulb pump, airflow becomes restricted and cushion stability, along with the ability to position to the prescribed seating posture is achieved. Each cell is effectively “LOCKED” into place. This unique “cell lock” technology allows you to address multiple seating challenges in an effective and efficient way.

Standard Contour Cushion

The Standard Contour Cushion (4" cell height) is Star Cushion Products premier cushion for the prevention and treatment of pressure ulcers. Flexible, interconnected 4" air cells allow for deep immersion into the cushion providing constant low pressures across the contact seating surface. In addition, the cells move and shape to the contact area providing a low shear, low friction surface ideal for the prevention and healing of up to and including a Stage IV pressure ulcer.



Product Features:

- Air adjustable 4" cells offer maximum immersion for the ultimate in pressure reduction, and tissue protection and treatment.
- Available in dual and four valve/compartments to effectively address seating challenges such as pelvic obliquities and sacral sitting, to enhance posture control and to improved overall function.
- Light weight, Approximately 3 pounds/1.36 kilograms
- Available in Pediatric, Adult, and Bariatric Sizes
- 24 Month Manufacturers Limited Warranty
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene rubber.

Mid contour Cushion




The Mid Contour Cushion (3" cell height) is recommended for individuals who require a cell height that provides a high level of tissue protection and treatment while offering greater ease in transfers in and out of the cushion. The Mid Contour cell height has proven to be ideal for the pediatric client, providing deep immersion for high levels of tissue protection and treatment along with a lower cell height for the pediatric frame.

Product Features:

- 3" high cell is recommended for those individuals who require deep immersion for tissue protection and treatment, but prefer the lower cell height for basic daily activities. Ideal for the pediatric client.
- Available in dual and four valve/compartments to effectively address seating challenges such as pelvic obliquities and sacral sitting, to enhance posture control and to improved overall function.
- Light Weight, Approximately 3 pounds/1.36 kilograms
- 24 Month Manufacturers limited Warranty
- Available in Pediatric, Adult and Bariatric Sizes
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene rubber.


Short Contour Cushion

The Short Contour Cushion (2" cell height) is recommended for the active individual that can perform independent weight shifts, or for those individuals who are at low risk for pressure ulcers. The Short Contour is the perfect cushion for sports enthusiasts as well, providing the ideal seating surface for your everyday chair, sports chair or other sport equipment.



Product Features:

- 2" cell cushion ideal for the individual who has an active lifestyle, is "at low risk" for, or has a history of pressure ulcers.
- Available in dual and four valve/compartments to effectively address seating challenges such as pelvic obliquities and sacral sitting, to enhance posture control and to improved overall function.
- Light Weight, Approximately 2.5 pounds/1.13 Kilograms
- 24 Month Manufacturers limited Warranty
- Available in Pediatric, Adult and Bariatric Sizes
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene rubber.



CXR Cushion


The CXR™ Cushion is designed to minimize coccyx pressures through a combination of 2" Cells under the coccyx area, and 4" cells throughout the remaining cushion surface. The dual valve configuration allows you to adjust both the 4" and 2" cells separately achieving individualized immersion, and ultimately constant low sitting pressures. The 2" cells allow for immersion while also providing added relief and protection under the coccyx area. 16 standard configurations address specific individualized seating needs. Ideal for the individual who has, or is at risk for the development of pressure ulcers.

Product Features:

- The CXR is designed for those individuals who required added protection and treatment under the coccyx. 16 standard configurations available to better meet the specific seating needs of the individual.
- Standard as a dual valve with both the 2" cells and the 4" cells being valved separately.
- Ideal for the minimization of coccyx pressures.
- Light Weight, approximately 3 pounds, 1.36 kilograms.
- 24 Month Manufacturers Limited Warranty.
- Available in Pediatric, Adult and Bariatric Sizes.
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene rubber.

STAR[★]MATT™

NON POWERED ADJUSTABLE ZONED PRESSURE REDUCING MATTRESS OVERLAY



The StarMatt Overlay is a four section, therapeutic, air adjustable mattress overlay featuring 180 air cells per section or 720 air cells in a total StarMatt system. As the individual is immersed into the StarMatt, they benefit from a high resolution, low pressure, low shear support surface that conforms to the unique contours of the body, providing the ideal environment for the healing of up to and including a Stage IV pressure ulcer. The StarMatt, the “**GOLD**” Standard in therapeutic mattress overlays.

A vapor permeable, quilted top cover is included with each StarMatt System. The cover encompasses not only the StarMatt but the hospital or bed mattress as well, protecting the bed mattress from fluids and helping to prevent sliding of the StarMatt during Transfers in and out of bed. Machine washable.



A heavy duty carrying bag is included with each StarMatt System providing a convenient way to transport the StarMatt from individual care settings. The carrying bag easily holds all four StarMatt sections, cover and hand pump.

Product Features:

- Four section mattress system fitting a standard size hospital bed.
- Each section is adjustable to accommodate the unique body shapes and sizes of each individual.
- Each section is designed with 180 individual air cells, or 720 total air cells in a four section system, providing a high resolution, low shear, low friction, low pressure environment ideal for the healing of pressure ulcers, up to and including a Stage IV.
- Light weight, approximately 32 pounds or 14.5 kilograms per total system weight.
- Non powered, so no pumps or powered components. Low maintenance
- 4” high cell for maximum immersion and pressure distribution.
- 18 month Manufacturers Limited Warranty
- Comes complete with four mattress sections, quilted top cover, carrying bag, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene rubber.

Function



Recliner Cushions, Commode Cushions and Custom Cushions



Product Features:

- Recliner cushions to fit standard recliner wheelchairs. Available in 2", 3" and 4" cell heights. Utilize a combination of cell heights for a customized fit on both seat and back cushion.
- Commode Seats and Shower chair cushions standard and available upon request.
- Custom Cushions, designed to fit your specific seating needs. Combine cell heights, design in cut-outs, add multiple valves or compartments.
- Quick turn-around times for Custom cushions. Custom cushions normally ship one week from receipt of order.
- Adjustable for deep immersion and constant low pressures. Conforms to unique shapes and sizes.
- Light Weight
- 24 Month Limited Manufacturers Warranty.
- Custom cushions and recliners come complete with fitted cover, hand pump and repair kit. Commode seat and shower chair cushions do not include a fitted cover.
- Designed of high quality, flame resistant, latex free neoprene rubber.

Item Numbers and Sizing Chart

* the "-1" following the item # indicates number of valves/compartments provided per cushion. Single, dual and four valve/compartment configurations available as standard items upon request and at no additional charge.

STANDARD CONTOUR™	MID CONTOUR™	SHORT CONTOUR™	THE CXR™	STARLOCK™	RECOMMENDED CHAIR SIZE
SC 1515-1	MC 1515-1	SCS 1515-1	CXR 1515-2	SL 1515-1	16" x16"
SC 1517-1	MC 1517-1	SCS 1517-1	CXR 1517-2	SL 1517-1	16" x18"
SC 1519-1	MC 1519-1	SCS 1519-1	CXR 1519-2	SL 1519-1	16" x20"
SC 1521-1	MC 1521-1	SCS 1521-1	CXR 1521-2	SL 1521-1	16" x22"
SC 1715-1	MC 1715-1	SCS 1715-1	CXR 1715-2	SL 1715-1	18" x16"
SC 1717-1	MC 1717-1	SCS 1717-1	CXR 1717-2	SL 1717-1	18" x18"
SC 1719-1	MC 1719-1	SCS 1719-1	CXR 1719-2	SL 1719-1	18" x20"
SC 1721-1	MC 1721-1	SCS 1721-1	CXR 1721-2	SL 1721-1	18" x22"
SC 1915-1	MC 1915-1	SCS 1915-1	CXR 1915-2	SL 1915-1	20" x16"
SC 1917-1	MC 1917-1	SCS 1917-1	CXR 1917-2	SL 1917-1	20" x18"
SC 1919-1	MC 1919-1	SCS 1919-1	CXR 1919-2	SL 1919-1	20" x20"
SC 1921-1	MC 1921-1	SCS 1921-1	CXR 1921-2	SL 1921-1	20" x22"
SC 2115-1	MC 2115-1	SCS 2115-1	CXR 2115-2	SL 2115-1	22" x16"
SC 2117-1	MC 2117-1	SCS 2117-1	CXR 2117-2	SL 2117-1	22" x18"
SC 2119-1	MC 2119-1	SCS 2119-1	CXR 2119-2	SL 2119-1	22" x20"
SC 2121-1	MC 2121-1	SCS 2121-1	CXR 2121-2	SL 2121-1	22" x22"



Exclusive Canadian Distributor

3223 Orlando Drive, Mississauga, ON, L4V 1C5, 905-671-1661, 888-737-4011, fax 866-670-3377 www.future-mobility.com
 Future Mobility Healthcare Inc. Orion and Prism are trademarks of Future Mobility Healthcare Inc.