


PRISM Seating & Positioning



The Super Cushion for Kids is a Junior size cushion created with three layers of polyurethane, high density and memory foam provides children with excellent pressure relief.

The anatomically designed firm foam base with its aggressive leg troughs and abductor contouring provides superior pelvic stability and leg positioning.

A central core of pressure relieving visco-elastic insert under the IT's allows for excellent immersion for bony prominences.



SPOTLIGHT ON

- High density layered foam for a consistent, long lasting product
- Standard with water resistant outer cover made of black meshtex urethane
- Waterfall front edge for foot propelling comfort
- Radius cut to match shape of back upholstery
- Pressure relieving
- Optional water resistant inner cover

WARRANTY

Two (2) years limited on Cushion and 90 days on cover.

PRODUCT WEIGHT

1.50 lbs

USER WEIGHT

150 lbs

CUSHION THICKNESS

A—Height at adductor = 3 3/4"

B—Height at abductor = 4"

C—Height at rear IT = 2 1/2"

D—Leg trough = 2 1/2"

E— Height at rear trochanter = 3 1/2"

Super Cushion for Kids

Stock Sizes—same day shipments

10"W x 10"D	KSC1010	16"W x 14"D	KSC1614
12"W x 12"D	KSC1212	16"W x 16"D	KSC1616
14"W x 14"D	KSC1414		
14"W x 16"D	KSC1416		

Non-stocked Sizes—shipments within 3 business days

10"W x 12"D	KSC1012	14"W x 10"D	KSC1410
10"W x 14"D	KSC1014	14"W x 12"D	KSC1412
10"W x 16"D	KSC1016	16"W x 10"D	KSC1610
12"W x 10"D	KSC1210	16"W x 12"D	KSC1612
12"W x 14"D	KSC1214		
12"W x 16"D	KSC1216		

