

**SUPER CHAIR FOR KIDS BY ORION**

**Order Form/Price List**

(effective March 31,2008)

**Account Information:**

**Date:** \_\_\_\_\_

Account #: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_  
 Province: \_\_\_\_\_  
 Postal: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 Comment: \_\_\_\_\_

PO# \_\_\_\_\_  
 Ship To: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_  
 Province: \_\_\_\_\_  
 Postal: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_

**1 FRAME TYPE**

- Tilting Frame w/Swingaway Front End
- 10"-13" Frame (crosst tubing)      SUPERCHAIR      \$3,250      CT10      No Charge
- 13"-16" Frame (crosst tubing)      CT13      No Charge

**FRAME FINISH**

- Black      FBK      No Charge
- Blue      FBE      No Charge
- Burgundy      FBY      No Charge
- Orange      FOE      No Charge
- Yellow      FYW      No Charge
- Red      FRD      No Charge
- Silver      FSR      No Charge

**SEAT WIDTH**

- 10" Seat Width      SW10      No Charge
- 11" Seat Width      SW11      No Charge
- 12" Seat Width      SW12      No Charge
- 13" Seat Width      SW13      No Charge
- 14" Seat Width      SW14      No Charge
- 15" Seat Width      SW15      No Charge
- 16" Seat Width      SW16      No Charge

**SEAT DEPTH**

- 11" Seat Depth      SD11      No Charge
- 12" Seat Depth      SD12      No Charge
- 13" Seat Depth      SD13      No Charge
- 14" Seat Depth      SD14      No Charge
- 15" Seat Depth      SD15      No Charge
- 16" Seat Depth      SD16      No Charge

**SEAT-TO-FLOOR HEIGHT - front - 20", 22" Rear Wheels**

- 2  17.5" (20" rear & 4" caster)      SFK01      No Charge
- 3  18" (20" rear & 4" caster)      SFK02      No Charge
- 18.5"(22" rear & 5" caster)      SFK03      No Charge
- 19" (22" rear & 7" caster)      SFK04      No Charge
- 19.5"(22" rear & 7" caster)      SFK05      No Charge
- 19.75"(22" rear & 8" caster)      SFK06      No Charge

**SEAT-TO-FLOOR HEIGHT - front - 12" Rear Wheels**

- 4  18.5" (12" rear & 5" caster)      SFK07      No Charge
- 19" (12" rear & 7" caster)      SFK08      No Charge
- 19.5" (12" rear & 7" caster)      SFK09      No Charge
- 20" (12" rear & 8" caster)      SFK10      No Charge

**Basic Chair Includes:**

Rigid Tilting Frame w/Swingaway Front End  
 Round Heavyduty Crossbrace Tubing  
 Seat Width (10"-16")  
 Seat Depth (11"-16")  
 Seat to Floor Height Front (17.50"-20")  
 Seat to Floor Height Rear (15.75"-18.5")  
 Rear Wheel Options (12", 20" or 22" )  
 Handrim Options (Stainless Steel, Plastic Coated)  
 Caster Options (4", 5", 7" or 8")  
 Adjustable Height Full Length Arm T-Style  
 Super Chair Kids Cushion  
 Super Chair Kids Back  
 Angle Adjustable Push Canes  
 Aluminum Flip-Up Footplate  
 Heel Loops  
 Wheel Locks (Push or Pull to Lock)  
 Side Storage Pocket  
 Padded Headrest (Kids or Small)  
 Choice of Seven Colours (Black, Blue, Burgundy, Orange, Yellow, Red or Silver)

**FOOTNOTES:**

1. Seat angle is set at 3 degrees. See service manual for details
2. Can not use 8" caster with 20" rear wheel and large fork
3. Can not use 4" caster with 22" rear wheel and small fork
4. Can not use 4" , 5" or 8" caster with 12" rear wheel and small fork
5. 12" Rear Wheel does not come with handrims
6. Super Chair Kids Back Stock 10"W, 12"W, 14"W, 16"W
7. Super Chair Kids Cushion Stock 10"W, 12"W, 14"W, 16"W
8. Low Cane 20 3/4", High Cane 24 3/4"
9. Super Chair Kids Aluminum Footplate available for chairs 10"W to 16"W
10. Front Seat Pan available in two depths 11" to 14" or 13" to 16"
11. Rear Seat Pan available in two depths 11" to 14" or 13" to 16"
12. Composite Flip-Up Footplates available for chairs 14" to 16" W
13. Aluminum Flip-Up Footplates available for chairs 12" to 16" W
14. Angle Adj. Flip-Up Footplates available for chairs 15" to 16" W
15. Kids Headrest Foam Pad available for chairs 10"W to 14"W
16. Small Headrest Foam Pad available for chairs 15"W & 16"W

**Future Mobility Healthcare Inc. reserves the right to change prices and specifications without notice.**

Note: All specifications and dimensions are approximate  
 Prism and Orion Products are trademarks of Future Mobility Healthcare Inc.

**WHEEL OPTIONS**

5	<input type="checkbox"/> 12" Composite Urethane Tires	WCK12	No Charge
	<input type="checkbox"/> 20" Composite Urethane Tires	WCK20	No Charge
	<input type="checkbox"/> 22" Composite Urethane Tires	WCK22	No Charge

**HANDRIM OPTIONS**

<input type="checkbox"/> 20" Stainless Steel	HSK20	Standard
<input type="checkbox"/> 22" Stainless Steel	HSK22	Standard
<input type="checkbox"/> 20" Plastic Coated	HPK20	\$160
<input type="checkbox"/> 22" Plastic Coated	HPK22	\$160
<input type="checkbox"/> Omit Handrims	OHR	less net \$15

**CASTER OPTIONS**

<input type="checkbox"/> 4" x 1" Urethane Tire	UTK4	No Charge
<input type="checkbox"/> 5" x 1.25" Urethane Tire	UTK5	No Charge
<input type="checkbox"/> 7" x 1.25" Urethane Tire	UTK7	No Charge
<input type="checkbox"/> 8" x 1.25" Urethane Tire	UTK8	No Charge

**ARM TYPE**

<input type="checkbox"/> Full Length Adj. Height (vinyl pad) T-Style	AFLK	\$125
<input type="checkbox"/> Desk Length Adj. Height (vinyl pad) T-Style	ADLK	\$125
<input type="checkbox"/> Curved Gel Adj. Height T-Style	ACGK	\$125
<input type="checkbox"/> Curved Foam Adj. Height T-Style	ACFK	\$125
<input type="checkbox"/> Cantilever Arm Full Length Adj. Angle (vinyl pad)	ACFLK	\$225
<input type="checkbox"/> Cantilever Arm Desk Length Adj. Angle (vinyl pad)	ACDLK	\$225
<input type="checkbox"/> Lower Arm Rest w/guard	ALOK	\$125

**8 BACK TYPE**

<input type="checkbox"/> Angle Adj 10 deg Bend Cane w/Push bar (low)	BLPK	\$150
<input type="checkbox"/> Angle Adj 10 deg Bend Cane w/Push bar (high)	BHPK	\$150
<input type="checkbox"/> Angle Adj 10 deg Bend Cane w/o Push bar (low)	BLNPK	\$150
<input type="checkbox"/> Angle Adj 10 deg Bend Cane w/o Push bar (high)	BHNPK	\$150

**BACK OPTIONS**

6	<input type="checkbox"/> Super Chair Kids Back Stock	BSK	Standard
	<input type="checkbox"/> Super Chair Kids Back Non Stock	BNK	Standard
	<input type="checkbox"/> Super Chair Kids Back Custom	BCMK	\$825
	<input type="checkbox"/> Slip on Padded Back Upholstery	SBUK	Standard
	<input type="checkbox"/> Prism Basic Back on chair Custom	BBCUS	\$825
	<input type="checkbox"/> Prism Ultra Back on chair Custom	UBCUS	\$1,010
	<input type="checkbox"/> Prism Air Back on chair Custom	ABCUS	\$1,075
	<input type="checkbox"/> Prism TrueFitt Back on chair Custom	TFBCUS	\$1,025
	<input type="checkbox"/> Omit Back Upholstery	OB	less net \$30
	<input type="checkbox"/> Specify Back Size W____x H____		

**SEAT OPTIONS**

10	<input type="checkbox"/> Front Seat Pan	FSP	Standard
11	<input type="checkbox"/> Rear Seat Pan	RSP	Standard
7	<input type="checkbox"/> Super Chair Kids Cushion Stock	CSK	Standard
	<input type="checkbox"/> Super Chair Kids Cushion Non Stock	CNK	Standard
	<input type="checkbox"/> Super Chair Kids Cushion Custom	CCMK	\$685
	<input type="checkbox"/> Prism Ideal Cushion - Custom	ICCM	\$275
	<input type="checkbox"/> Prism Elite Cushion - Custom	ECCUS	\$425
	<input type="checkbox"/> Prism Supreme Cushion - Custom	SCCUS	\$685
	<input type="checkbox"/> Prism Supreme Air Cushion - Custom	SACCUS	\$785
	<input type="checkbox"/> Prism Supreme II Plus Cushion - Custom	SCPCUS	\$685
	<input type="checkbox"/> Omit Cushion	OC	less net \$10

**FRONT RIGGING**

<input type="checkbox"/> 70 Degree Swingaway Upgraded	A70K	No Charge
<input type="checkbox"/> 80 Degree Swingaway Upgraded	A80K	\$215
<input type="checkbox"/> 90 Degree Swingaway Upgraded	A90K	\$230
<input type="checkbox"/> Elevating Legrest with Swingaway Upgraded	ELRK	\$550
<input type="checkbox"/> Omit Front Rigging	OFR	less net \$20

**FOOTPLATE TYPE**

9	<input type="checkbox"/>	Super Chair Kids Aluminum Flip-Up Footplate	KAF	Standard
12	<input type="checkbox"/>	Composite Flip-Up Footplate	FPC	Standard
13	<input type="checkbox"/>	Aluminum Flip-Up Footplate	FPA	Standard
14	<input type="checkbox"/>	Adjustable Angle Flip-Up Footplate	FFAA	\$120

**FOOTREST/LEGREST ACCESSORIES**

<input type="checkbox"/>	Padded Calf Strap	PCS	\$70
<input type="checkbox"/>	Heel Loops	HL	Standard
<input type="checkbox"/>	Omit Heel Loops	OHL	less net \$7

**WHEEL LOCKS**

<input type="checkbox"/>	Push-to-Lock Wheel Lock	PSK	Standard
<input type="checkbox"/>	Pull-to-Lock Wheel Lock	PLK	No Charge
<input type="checkbox"/>	Wheel Lock Extension	WLK	No Charge

**HEADREST HARDWARE OPTIONS**

<input type="checkbox"/>	Fixed Headrest Hardware	OFH1	Standard
<input type="checkbox"/>	Single Multi-Axial Headrest Hardware	OFH2	\$158
<input type="checkbox"/>	Dual Multi-Axial Headrest hardware	OFH3	\$195
<input type="checkbox"/>	Fixed Headrest Hardware Ass'y (+ 4")	OFH4	\$20
<input type="checkbox"/>	Single Multi-Axial Headrest Hardware (+4")	OFH5	\$183
<input type="checkbox"/>	Dual Multi-Axial Headrest hardware (+ 4")	OFH6	\$220

**HEADREST FOAM PAD OPTIONS**


15	<input type="checkbox"/>	Kids Headrest Foam Pad	KHF	Standard
16	<input type="checkbox"/>	Small Headrest Foam Pad	OHP1	Standard
	<input type="checkbox"/>	Small Headrest with Wings Foam Pad	OHP2	\$90
	<input type="checkbox"/>	Large Headrest Foam Pad	OHP3	No Charge
	<input type="checkbox"/>	Large Headrest with Wings Foam Pad	OHP4	\$95
	<input type="checkbox"/>	Omit Headrest	OHT	less net \$10

**MISCELLANEOUS ACCESSORIES**


<input type="checkbox"/>	Anti-tippers	ATS	Standard
<input type="checkbox"/>	Auto-Seat Positioning Straps	ARB	\$52
<input type="checkbox"/>	Velcro Release Positioning Strap	VRB	\$48
<input type="checkbox"/>	4 Point Harness	FPH	\$95
<input type="checkbox"/>	Quick Release Axles	QRA	\$110
<input type="checkbox"/>	Front Transit Hook Assembly	FTHK	Standard
<input type="checkbox"/>	Rear Transit Hook Assembly	RTHK	Standard
<input type="checkbox"/>	Padded Lap Tray	PLT	\$285
<input type="checkbox"/>	Half Padded Lap Tray, Right	SHPLTL	\$145
<input type="checkbox"/>	Half Padded Lap Tray, Left	SHPLTR	\$145

**SEAT TO FLOOR HEIGHT**


STF	REAR	CASTER	FORK	FORK HOLE	AXLE PLATE	STF
FRONT	WHEEL	SIZE	SIZE		HOLE	REAR
(inches)	(inches)	(inches)				(inches)
18.50	12	4	LG	C	A	17.25
18.50	12	5	LG	B	A	17.25
19.00	12	7	SML	C	A	17.50
19.50	12	7	LG	B	B	18.25
20.00	12	8	LG	B	B	18.50
17.50	20	4	SML	C	A	15.75
18.00	20	5	SML	C	A	15.75
18.50	22	5	LG	B	A	16.75
19.00	20	7	SML	C	B	17.00
19.50	22	7	LG	B	B	18.00
20.00	22	7	LG	C	B	18.00




**LONG FORK**



**SHORT FORK**



**12" WHEEL AXLE PLATE**



**20", 22" WHEEL AXLE PLATE**